

# Combination agreement Kombinationsavtal

Wheelchairs/Rullstolar



1A15211.13 24-07-01

**en** The list of combination agreements is only an overview. For more detailed information regarding the agreements, please contact our customer service. The maximum user weight for the combination is always the same as for the product in the combination that has the lowest maximum user weight.

**sv** Listan över kombinationsavtal är endast en översikt. För mer ingående information angående avtalen, kontakta kundservice. Den maximala brukarvikten för kombinationen är alltid densamma som för den produkt i kombinationen som har den lägsta maximala brukarvikten.

Manufacturer/Tillverkare	Product/Produkt	Etac Product/Produkt
Alber	e-motion M15, M25	Cross 6, Elite, Cross, Cross XL, Cross 5, Cross 5 XL, M100, Twin, Prio, Crissy Swing-Away, Crissy Active
Alber	e-fix E25/E26	Cross 6, Cross, Cross XL, Cross 5, Cross 5 XL Transit, Twin, M100
Alber	Quix Q10	Cross, Cross XL, Transit, Twin
Alber	Scalamobil S30 equipped with Scalaport X7	Act, Elite, Next, Next Comfort, M100
Alber	Scalamobil S35/S38	Cross 6, Cross, Cross XL, Cross 5, Cross 5 XL, Transit, Twin, M100, Prio, Crissy Swing-Away, Crissy Active
Alber	Viamobil eco V14	Cross, Cross XL, Cross 5, Cross 5 XL, Twin, M100, Crissy Swing-Away, Crissy Active
Alber	Viamobil V15	Cross, Cross XL, Transit, Twin, M100
Alber	Viamobil V25	Cross 6, Cross, Cross XL, Cross 5, Cross 5 XL, M100, Twin, Transit, Prio, Crissy Swing-Away, Crissy Active
Alber	e-fix E35/E36	Cross 6, Cross, Cross XL, Cross 5, Cross 5 XL, M100, Twin, Transit, Prio, Crissy Swing-Away, Crissy Active
Alber	e-pilot P15	Cross 6, Crissy Swing Away, Crissy Active
Alber	SMOOV 010	Cross 6, Cross 5, Cross 5 XL, Crissy Swing Away, Crissy Active
Alu Rehab A/S	Netti Contour, UNO seat, Netti Sit, UNO Back, Netti Smart, Netti Stabil, Netti Super Stabil, Netti Kyfotic Back, Netti Nordic Back.	Cross, Cross XL, Cross 5, Cross 5 XL, Act, Elite, M100, M100 TR, Transit III, Next, Next Comfort, Sting.
Anatomic Sitt	Proserie Headrest	Cross 6, Next, Cross, Cross XL, Twin, Transit III, Sting, Act, Elite, Crissy Swing Away, Crissy Active, Cross 5, Cross 5 XL, Prio
Armbike/Jonny Borg	Armbike	Act, Elite
Bodypoint, Inc	Bodypoint 2-pt and 4-pt Hip Belts: NB200-series, HB100-series, HB200-series, HB400-series, EB200-series, LB200-series, LB400-series, LH600-series, SH100-series, SH200-series, SH300-series, SP110-series, SP111-series	Cross 6, Cross, Cross XL, Cross 5, Cross 5 XL, Twin, Transit, M100, M100 TR, Prio
Decon AB	Easy Wheel	Prio, Cross, Twin, Act, Elite, Xact
Decon AB	E-Move, E-Drive	Prio, Cross 6, Cross 5, Cross 5XL, Cross 5 Active, Xact with 6" castors
Decon AB	Adventus, Infinity	Xact

# Combination agreement Kombinationsavtal

Wheelchairs/Rullstolar



1A15211.13 24-07-01

Manufacturer/Tillverkare	Product/Produkt	Etac Product/Produkt
Enkatsu Solutions AB	Enkatsu Climber	Cross, Twin, Transit, Ono Rollator, Clean, Mobil
G.Niklassons	Maxgrip Rear wheel	Next, Cross XL, Cross, Twin, Sting, Act, Elite
HEA Medical AB	HEA attachment for chest slings	Cross 6, Cross 5, Cross 5 XL, Cross, Cross XL, M100 with height adjustable push handles, Transit with height adjustable push handles
Max mobility	Smart Drive MX2+	Cross 5, Cross 5 XL, Cross, Cross XL
Medema	M-Drive power assist add on	Cross, Cross XL, Cross VB & Cross XL VB
Mercado Medic	Real 9000, Real 9000 Plus, Real 6100, Real 6100 Plus	Amputee leg support, Leg support- wide -narrow and -narrow short, Foot supports
Otto Bock Scandinavian AB	Head and Neck supports, Foot Buckets and Calf Pads, Firm Seating and Back systems and Cushions	Next, Cross XL, Cross, Twin, Transit III
Performance Health	Backrest V-Trak Products LTD	Cross, Cross XL, Next, Next Comfort, Twin, Transit
R82	Head and Necks support	Cross, Cross XL, Cross 5, Next, Next Comfort, Elite, Sting, M100, M100 TR
Sunrise Medical AB	Backrest Jay 2, standard, deep & Back Jay 3	Cross 6, Cross 5, Cross 5 XL, Prio, Cross, Cross XL, Next, Next Comfort, Twin, Transit
Sunrise Medical AS	PowerTec F16	Cross
Sunrise Medical Europe	Quickie Wheeldrive, Empulse R20	Cross 5, Cross 5 XL, Cross 6, Prio, Crissy
Särnmark Hjälpmedel AB	Hemi-armrest Positioning seat & Active badge	Etac´s wheelchair
Tzora Active Systems	Combi Power, Combi Power Comfort	Cross, Cross XL, Transit III
UB-Let	EasyGo+	Cross XL, Twin XL
	Seat cushions, CE-marked acc. to MDR and intended for wheelchair use	Etac manual wheelchairs
	Belts and soft trunk supports, CE-marked acc. to MDR and intended for wheelchair use	Etac manual wheelchairs