

Combination agreement

Wheelchairs

KLAXON
by Etac

1A15211.15 24-11-11

en The list of combination agreements is only an overview. For more detailed information regarding the agreements, please contact our customer service.
The maximum user weight for the combination is always the same as for the product in the combination that has the lowest maximum user weight.

Manufacturer	Product	Klaxon Product
Aria Wheels	2.0A1, 1.0Mg, ULTRA, KID	KLICK Electric (Standard / Tetra), KLICK Power (Standard / Tetra), KLICK Race (Standard / Tetra), KLICK Mini (Standard / Carbon), KLICK Monster, KLICK Hybrid
Invacare	Küschall Compact version 2.0 (lateral fixation), Küschall K-Series version 2.0 (central fixation), Küschall The KSL version 2.0 (lateral fixation and trapezoidal clamps), Küschall Champion version 2.0 (lateral fixation and trapezoidal clamps)	KLICK Power (Standard / Tetra), KLICK Race (Standard / Tetra), KLICK Hybrid
Merya	NANO Model 1.155, NANO X Model 1.156, Smart F Model 2.360, Smart S Model 2.370, Avanti Pro Model 1.735, Avanti Model 1.736, Hurricane Model 1.880	KLICK Electric (Standard / Tetra), KLICK Power (Standard / Tetra), KLICK Race (Standard / Tetra), KLICK Mini (Standard / Carbon), KLICK Monster, KLICK Hybrid, KLICK Power Limited Edition (Standard / Tetra), KLICK Power (2021), KLICK Race (2021), KLICK Monster (2021)
Motion Composites	MOVE, A6, A7, C2, VELOCE, APEX Aluminum, APEX Carbon, Hello Kids	KLICK Electric (Standard / Tetra), KLICK Power (Standard / Tetra), KLICK Race (Standard / Tetra), KLICK Mini (Standard / Carbon), KLICK Monster, KLICK Hybrid
Off Carr	EOS, QUASAR, VENUS, FENICE, ALHENA, DIVA, ALTHEA, QPX, ZODIAC, IDRA 2.0, HALLEY, VEGA	KLICK Electric (Standard / Tetra), KLICK Power (Standard / Tetra), KLICK Race (Standard / Tetra), KLICK Mini (Standard / Carbon), KLICK Monster, KLICK Hybrid
Otto Bock	Avantgarde 4, Zenit, Ventus, Motus, Invader, Start 4 M1, Start 4 M2, Start 4M2S	KLICK Electric (Standard / Tetra), KLICK Power (Standard / Tetra), KLICK Race (Standard / Tetra), KLICK Mini (Standard / Carbon), KLICK Monster, KLICK Hybrid
Rehateam Progeo	Exelle Vario, Yoga, Ego, Ego Custom, Ego Carbon, Joker Standrad, Joker Full Carbon, Joker Carbon Pack, Joker Evolution, Joker R2, Carbomax, Duke, Noir 1.0, Noir 2.0, Tekna Advance, Tekna Advance Carbon	KLICK Electric (Standard / Tetra), KLICK Power (Standard / Tetra), KLICK Race (Standard / Tetra), KLICK Mini (Standard / Carbon), KLICK Monster, KLICK Hybrid
RGK Wheelchairs	Tiga, Tiga Sub4, Tiga FX, Tiga F2, Octane, Octane Sub4, Octane FX, Hilite, Hilite XTR, Maxlite	KLICK Electric (Standard / Tetra), KLICK Power (Standard / Tetra), KLICK Race (Standard / Tetra), KLICK Monster, KLICK Hybrid, Klaxon Twist
Sunrise Medical	Helium, Nitrum, Argon/Argon ² , Neon, Neon ² , Life, Life R, Xenon ² , Easy 200/300, Helium	KLICK Electric (Standard / Tetra), KLICK Power (Standard / Tetra), KLICK Race (Standard / Tetra), KLICK Mini (Standard / Carbon), KLICK Monster, KLICK Hybrid, KLICK Power Limited Edition (Standard / Tetra), KLICK Power (2021), KLICK Race (2021), KLICK Monster (2021)
Vassilli	Folding Wheelchairs: 17.60N, 17.58N, 17.57N, 17.70N Rigid Wheelchairs: 18.64M and 12.70L	KLICK Electric (Standard / Tetra), KLICK Power (Standard / Tetra), KLICK Race (Standard / Tetra), KLICK Mini (Standard / Carbon), KLICK Monster, KLICK Hybrid

Combination agreement

Wheelchairs

KLAXON
by Etac

1A15211.15 24-11-11

Manufacturer	Product	Klaxon Product
Vermeiren Group	Jazz S50, V300, D200, D100, D200P, Eclips X4, Eclips+, V100, D20030°, V500, V500XR, V300XR, Trigo S, Trigo T	KLICK Electric (Standard / Tetra), KLICK Power (Standard / Tetra), KLICK Race (Standard / Tetra), KLICK Mini (Standard / Carbon), KLICK Monster, KLICK Hybrid, KLICK Power Limited Edition (Standard / Tetra), KLICK Power (2021), KLICK Race (2021), KLICK Monster (2021)
Wolturnus	W5, W5 A, W5 S, W5 K, W5 SL, W5 D, W5 XXL Tukan A, Tukan B, Tukan C, Tukan D, Tukan Tilt Hawk Dalton, Dalton Low, Dalton F, Dalton Light Merlin, Merlin K, Merlin D	KLICK Electric (Standard / Tetra), KLICK Power (Standard / Tetra), KLICK Race (Standard / Tetra), KLICK Mini (Standard / Carbon), KLICK Monster, KLICK Hybrid