

# Convoid Carrot 3 Booster Seat - Child Restraint System

## Canadian Price List / Order Form

Account #: \_\_\_\_\_ Purchase Order #: \_\_\_\_\_ Date: \_\_\_\_\_

Contact: \_\_\_\_\_ ATP/RTS#: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

### Carrot 3 Booster Seat Standard Features:

- |  |  |
|--|--|
| <ul style="list-style-type: none"> <li>• Carrot 3 Backrest</li> <li>• Recline Adjustment (90° - 180°)</li> <li>• Head Support: 8.5" x 10.5" (H x L)</li> <li>• 4 Point Positioning Harness</li> <li>• H-Harness Attaching Bracket</li> </ul> | <ul style="list-style-type: none"> <li>• Booster Seat Base</li> <li>• Lower Anchors Connector System</li> <li>• 2" Backrest Extension</li> <li>• Shoulder Protector Wings</li> <li>• Tether Strap (Black)</li> </ul> |
|--|--|

#### Standard Models:

- |   |                |
|---|----------------|
| <input type="checkbox"/> Carrot 3 Booster Seat CRT3BLACK - CCRT3BS-10 | <b>Price</b>   |
| <input type="checkbox"/> Carrot 3 Booster Seat CRT3BLUE - CCRT3BS-20  | <b>\$4,284</b> |
| <input type="checkbox"/> Carrot 3 Booster Seat CRT3RED - CCRT3BS-30   | <b>\$4,284</b> |

All prices in CAD

#### Carrot Booster Seat Recommended Item

- |   |              |
|---|--------------|
| <input type="checkbox"/> Top Securement Strap*        | <b>Price</b> |
| <input type="checkbox"/> Lower Securement Strap*      | <b>\$79</b>  |
| <b>*For vehicles not equipped with anchor points.</b> | <b>\$99</b>  |

Black	Blue	Red
CCT3HSF-10	CCT3HSF-10	CCT3HSF-10
CCT3LFS-10	CCT3LFS-10	CCT3LFS-10

#### Accessories

- |   |              |
|---|--------------|
| <input type="checkbox"/> Booster Support Tray (Crash Tested)                      | <b>Price</b> |
| <input type="checkbox"/> 3D Body Balance Pads L/R Set                             | <b>\$438</b> |
| <input type="checkbox"/> Seat Angle Adjustment Wedge (Height: 4.25" / Angle: 20°) | <b>\$227</b> |
| <input type="checkbox"/> Backrest Extension - 4" (Extends back height to 32.8")   | <b>\$171</b> |
| <input type="checkbox"/> Large Crotch Belt  | <b>\$137</b> |
|   | <b>\$359</b> |

Black	Blue	Red
CCRT3BT-19	CCRT3BT-29	CCRT3BT-39
CCT3BBP-00	CCT3BBP-00	CCT3BBP-00
CCT3SA0-10	CCT3SA0-20	CCT3SA0-30
CCT3BE4-10	CCT3BE4-20	CCT3BE4-30
CCT3CBX-00	CCT3CBX-00	CCT3CBX-00

Standard Shipping Rate: Ground \$71.95 • For expedited shipping rates, please call 1-844-876-6245 or email [convoidsales.us@etac.com](mailto:convoidsales.us@etac.com)

Prices in CAD & subject to change, ex works Torrance, California.


Effective: 2/1/2025

Convoid  
2830 California Street  
Torrance, CA 90503  
1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)

OF1101-ENG-REV.01



**CLIENT MEASUREMENTS:**



**Carrot 3 Base Model Specifications:**

Seat Width:	16.5"
Seat Depth:	15" - 16.5" *
Backrest Height Adjustment:	28" - 33"
Head Support Height:	8.5"
Backrest Height:	29" - 34"
Backrest Extension Height (Standard):	2"
Backrest Extension Height (Option):	4"
Shoulder Wing Height (Standard):	6"
Shoulder Wing Width Adjustment:	15.25" - 19.25"
Product Weight:	21 lbs.
Recline Adjustment:	90° - 180°
Product Overall Width:	20" - 22"

**User Measurement Guide:**

Minimum User Weight:	79 lbs.
Maximum User Weight:	165 lbs.
Minimum User Height:	55"
Maximum User Height:	70"

\*Factory Setting