

# Convaid Coaster - WM0002143

## Canadian Price List / Order Form

Account #: \_\_\_\_\_ Purchase Order #: \_\_\_\_\_ Date: \_\_\_\_\_  
 Contact: \_\_\_\_\_ ATP/RTS: \_\_\_\_\_ Mark For: \_\_\_\_\_

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

### Coaster Standard Features:

- Three fixed tilt setting options 10°, 20°, or 30° (Adjustable with Tools)
- Four Fixed Recline Setting Options (Adjustable with Tools)
- Tension Adjustable Back Straps
- 2-Point Pelvic Belt
- Seat Depth Growth of 2"
- Standard Convaid Upholstery
- Zippered Pockets in Seat & Back Upholstery
- Removable, Height Adjustable Swing-Away Footplates
- One Piece Height Adjustable Push Handle

### Standard Models:

- Coaster CTR12 - C800121
- Coaster CTR14 - C800141
- Coaster CTR16 - C800161
- Coaster CTR18 - C800181

- 80° Fixed Front Rigging
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tires
- Quick Release Rolko Rear Wheels
- Single Foot Operated Wheel Locks
- Clamshell Compact Fold
- Closure Strap
- Saddle Bags
- Silver Powder Coated Frame

### Coaster Scout Additional Features:

- Height and Tension Adjustable Back Positioner
  - Height Adjustable, Flip-Up Armrest
  - Attendant Handbreak
  - 8" x 2" Front Ribbed, 12.5 x 2.5" Rear Pneumatic Tires
- Total for Scout Standard Features \$1,130**

**All Prices in CAD**

### Price

- \$2,775**
- \$2,775**
- \$2,775**
- \$2,775**

### Coaster Scout Models:

- Coaster Scout CTR12 - C800123
- Coaster Scout CTR14 - C800143
- Coaster Scout CTR16 - C800163
- Coaster Scout CTR18 - C800183

### Price

- \$3,904**
- \$3,904**
- \$3,904**
- \$3,904**

### Select Tilt and Recline Setting

Factory Default Setting if No Selection is Made: 10° Fixed Tilt and 105° Recline

#### Fixed Tilt set at 10°

Fixed Recline set at:  100 |  105 |  110 |  115

#### Fixed Tilt set at 20°

Fixed Recline set at:  90 |  95 |  100 |  105

#### Fixed Tilt set at 30°

Fixed Recline set at:  80 |  85 |  90 |  95



### Upholstery Options

#### Colors:

Panther Black

C903854

Royal Blue

C903850

Sassy Purple

C903851

Princess Pink

C903759

Forest Green

C903857

Cherry Red

C903855



#### Foam Kit:

Positioning Wedges (Foam Kit for Zippered Upholstery Pockets)

#### Lumbar Support:

Height and Tension Adjustable Back Positioner

#### Price

**\$141**

**\$537**

	CTR12	CTR14	CTR16	CTR18
<input type="checkbox"/>	C905567	C905568	C905569	C905570
<input type="checkbox"/>	C905538	C905618	C905539	C905619

**Seat Depth Extensions^**

- 3" Seat Depth Extension

**^The Coaster will ship with the +3 Seat depth extension that will only accomplish the +3 seat depth ranges (shown below).**

Price	CTR12	CTR14	CTR16	CTR18
\$141	C905540	C905540	C905638	C905638

**Future Seat Depth Growth Kit\***

- 3" Seat Depth Extension with Fixed Front Riggings (Growth Kit)^
- 3" Seat Depth Extension with Elevating Legrest (Growth Kit)^

**\*The Coaster will ship with standard seat depth and include a seat depth growth kit for future growth that will accomplish all seat depth ranges (shown below).**

**^Seat depth growth kit selection is dependent on the front rigging's selection.**

Price	CTR12	CTR14	CTR16	CTR18
\$141	C905639	C905639	C905640	C905640
\$141	C905641	C905642	C905643	C905644

**Seat Depth Ranges**

Seat Width	Standard 2"	Option +3"
12"	10-12	13-15
14"	11-13	14-16
16"	12-15	15-18
18"	14.5-17.5	17.5-20.5



**Fixed Front Riggings**  
Standard



**Elevating Legrest**  
Option

**Frame Options**

- Height Adjustable Flip-Up Armrests
- Caster Lock for 2" Wide Front Tires^
- Rear Anti-Tip Tubes (Cannot use with Curb Tipper)
- Curb Tipper (Cannot use with Rear Anti-Tip Tubes)

**^Manufacturer Installation Required.**

Price	CTR12	CTR14	CTR16	CTR18
\$344	C905611	C905611	C905611	C905611
\$141	C905526	C905526	C905526	C905526
\$174	C905527	C905527	C905527	C905527
\$148	C510737	C510737	C510737	C510737

**Transit Required Items\***

- Transit Anchor Brackets^
- Headrest Extension
- H-Harness with Padded Covers (Requires 3-Point Pelvic Belt)
- 3-Point Pelvic Belt
- Foot Positioners, Small 7"-10"\*
- Foot Positioners, Medium 8"-11"\*
- Foot Positioners, Large 9"-12"\*

**^Manufacturer Installation Required.**

**\*Crash Tested - WC19 and ISO 7176-19 compliant**

Price	CTR12	CTR14	CTR16	CTR18
\$302	C905531	C905531	C905531	C905531
\$195	C902863	C900484	C900138	C905571
\$148	C903988	C903988	C903988	C903988
\$73	C904017	C904017	C904017	C904017
\$236	C903415	C903415	C903415	C903415
\$236	C903416	C903416	C903416	C903416
\$236	C903417	C903417	C903417	C903417

**Transit Option**

- Q'Strain Transit Lap Belt

Price	CTR12	CTR14	CTR16	CTR18
\$279	C903094	C903094	C903094	C903094

**Tire Options**

- 7.5" x 2" Front Ribbed, 12.5" x 2" Rear Solid Knobby Tires
- 7.5" x 2" Front, 12.5" x 2" Rear Pneumatic Tires

Price	CTR12	CTR14	CTR16	CTR18
\$450	C905535	C905535	C905535	C905535
\$450	C905561	C905561	C905561	C905561

**Scout Option**

**Scout Standard Features:**

- Height and Tension Adjustable Back Positioner
- Height Adjustable, Flip-Up Armrest
- Attendant Handbreak
- 8" x 2" Front Ribbed, 12.5 x 2.5" Rear Pneumatic Tires

- 8" x 2" Front Ribbed, 12.5" x 2" Rear Solid Knobby Tires

Price	CTR12	CTR14	CTR16	CTR18
\$460	C905545	C905545	C905545	C905545
\$344	C905545	C905545	C905545	C905545
\$86	C905545	C905545	C905545	C905545
\$240	C905545	C905545	C905545	C905545
\$450	C905545	C905545	C905545	C905545

**Head Supports**

- Padded Headwing
- Occi Headwings
- Headrest Extension

Price	CTR12	CTR14	CTR16	CTR18
\$367	C905635	C905636	C905637	C905637
\$344	C905271	C905272	C905272	C905272
\$195	C902863	C900484	C900138	C905571

**Trunk Positioning**

- H-Harness with Padded Covers (Requires 3-Point Pelvic Belt)
- Full Torso Support Vest
- Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1")\*
- Side Panels

**\*Requires Side Panels**

Price	CTR12	CTR14	CTR16	CTR18
\$148	C903988	C903988	C903988	C903988
\$398	C905598	C905599	C905600	C905601
\$577	C904211	C904110	C904110	C904110
\$319	C905552	C905553	C905554	C905555

**Pelvic Positioning**

- 3-Point Pelvic Belt

Price	CTR12	CTR14	CTR16	CTR18
\$73	C904017	C904017	C904017	C904017

**Foot & Leg Positioning**

- Medial Thigh Support, Abductor
- Lateral Thigh Support, Adductor
- Calf Panel\*
- Elevating Leg Rest (90° - 180°)
- Padded Elevating Legrest Cover
- Individual Angle Adjustable Footplates (Elevating Leg Rest Required)
- Angle Adjustable Footplates
- One Piece Angle Adjustable Flip-Up Footplate (-15° - +15°)
- Fully Padded Footbox (One Piece Footplate Required)
- Footplate Securement Strap
- Foot Positioners, Small 7"-10"
- Foot Positioners, Medium 8"-11"
- Foot Positioners, Large 9"-12"

**\*Calf Panel cannot be used when footplate extension tube is rotated, and in the highest position.**

Price	CTR12	CTR14	CTR16	CTR18
\$244	C905576	C905577	C905578	C905579
\$244	C905580	C905581	C905582	C905583
\$93	C904580	C904580	C904580	C904580
\$377	C905541	C905542	C905543	C905544
\$141	C904619	C904619	C904619	C904619
\$240	C905624	C905625	C905626	C905627
\$240	C903445	C903268	C903268	C903269
\$167	C905557	C905558	C905559	
\$490	C904954	C904955	C904956	
\$47	C903584	C903584	C903584	C903584
\$236	C903415	C903415	C903415	C903415
\$236	C903416	C903416	C903416	C903416
\$236	C903417	C903417	C903417	C903417

**Medical Support Accessories**

- Medical Necessity Storage Basket
- IV Pole (Height Adjustable and Collapsible)<sup>^^ \*</sup>
  - Adj. Swing Away Hardware (Required for IV Pole)

**Price**  
**\$445**  
**\$288**  
**\$207**

CTR12	CTR14	CTR16	CTR18
C905587	C905588	C905588	C905589
C904566	C904566	C904566	C904566
C905612	C905612	C905612	C905612

**^^Cannot be used with canopy.**  
**\*Rear Anti-tip Tubes Required.**

**Accessories**

- Headrest Cover (Canopy)
- Extended Headrest Cover (Canopy)
- Under Seat Storage Basket
- Utility Bag
- Rain Cover (Canopy Required)
- Mosquito Net (Canopy Required)

**Price**  
**\$344**  
**\$377**  
**\$172**  
**\$64**  
**\$227**  
**\$227**

CTR12	CTR14	CTR16	CTR18
C905017	C905018	C905019	C905020
C905602	C905603	C905604	C905605
C905584	C905585	C905585	C905586
C904531	C904531	C904531	C904531
C905590	C905591	C905592	C905593
C905594	C905595	C905596	C905597

The Coaster transit option includes four manufacturer installed wheelchair transportation anchors. The Coaster transit model has been crash tested and conforms

*HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.*

**Standard Shipping Rate: \$71.95 • For expedited shipping rates, please call 1-844-876-6245 or email [convoidsales.us@etac.com](mailto:convoidsales.us@etac.com)**

Prices in CAD & subject to change, ex works Torrance, California.

Effective: 1/27/2025

Convoid  
 2830 California Street  
 Torrance, CA 90503  
 1 844-US-MOBILITY (844-876-6245)  
[www.etac.us.com](http://www.etac.us.com)

OF0201-ENG-REV.01

