

Convoid Safari Tilt - HCPCS E1232

Retail Price List / Order Form

Account #: _____

Purchase Order #: _____ Date: _____

Contact: _____ ATP/RTS: _____

Mark For: _____

Bill To:	Ship To:
Address:	Address:
City/State/Zip:	City/State/Zip:
Phone:	Phone:
Email:	

Safari Tilt Standard Features:

- 5° to 45° of Adjustable Tilt
- Silver Powder Coated Frame
- 7.5 x 2" Front, 11.5 x 2.5" Rear Solid Tires (SFT16)
- 6" Front, 11" Rear Solid Tires (SFT18)
- Quick Release Rear Wheels (SFT16)
- Foot Operated Wheel Locks
- 2-Point Pelvic Belt
- Heavy-Duty Back Upholstery
- One Piece, Non-Folding/Non-Elevating Footplate
- Support Cushion with Solid Seat Insert

Standard Models:

- | | |
|---|--------------------------------|
| <input type="checkbox"/> Safari SFT16 - C900526 | Price
\$3,111 |
| <input type="checkbox"/> Safari SFT18 - C902432 | Price
\$3,111 |

Heavy Duty Upholstery Colors

- | | |
|--|---------------------|
| <input type="checkbox"/> Panther Black | Price
N/C |
| <input type="checkbox"/> Royal Blue | Price
N/C |
| <input type="checkbox"/> Sassy Purple | Price
N/C |

SFT16	SFT18
C903854	C903854
C903850	C903850
C903851	C903851

Frame Options

- | | | |
|---|------------------------------|--------------|
| <input type="checkbox"/> Anatomic Back ¹ | Price
\$450 | HCPCS |
| <input type="checkbox"/> Recline Option ¹ | Price
\$620 | E1225 |
| <input type="checkbox"/> Height Adjustable Flip-Up Armrests | Price
\$307 | E0973 |
| <input type="checkbox"/> Attendant Hand Brakes ¹ (Knobby Tires Required) | Price
\$317 | |
| <input type="checkbox"/> Curb Tipper | Price
\$107 | |

SFT16	SFT18
C901319	C902476
C901321	
C901322	C902477
C904064	C904064
C904063	C904063

¹Manufacturer Installation Required.

Transit Required Items

- | | | |
|---|------------------------------|------------------------------|
| <input type="checkbox"/> Transit Anchor Brackets ^{1*} | Price
\$269 | HCPCS
K0108 |
| <input type="checkbox"/> Headrest Extension* | Price
\$141 | E0966 |
| <input type="checkbox"/> H-Harness with Padded Covers* (Requires 3-Point Pelvic Belt) | Price
\$107 | E0960 |
| <input type="checkbox"/> 3-Point Pelvic Belt* | Price
\$53 | E0978 |
| <input type="checkbox"/> Elevating Legrests with Swing-Away Footplates* | Price
\$363 | E0990 |
| Foot Positioners* - Select Size Below: | | |
| <input type="checkbox"/> Small 7" - 10" | Price
\$170 | K0108 |
| <input type="checkbox"/> Medium 8" - 11" (Default if size is not selected) | Price
\$170 | K0108 |
| <input type="checkbox"/> Large 9" - 12" | Price
\$170 | K0108 |

SFT16	SFT18
C900527	C902435
C901310	C902471
C903988	C903988
C904017	C904017
C901299	C902460
C903415	C903415
C903416	C903416
C903417	C903417

¹Manufacturer Installation Required.

*Crash Tested - WC19 and ISO 7176-19 Compliant

Transit Options

- | | | |
|---|------------------------------|------------------------------|
| <input type="checkbox"/> Q'Straint Transit Lap Belt | Price
\$202 | HCPCS
E0978 |
|---|------------------------------|------------------------------|

SFT16	SFT18
C903095	C903095

Head Supports

- | | | |
|---|------------------------------|------------------------------|
| <input type="checkbox"/> Padded Headwings | Price
\$265 | HCPCS
E0955 |
| <input type="checkbox"/> Occi Headwings | Price
\$249 | E0955 |
| <input type="checkbox"/> Headrest Extension | Price
\$141 | E0966 |

SFT16	SFT18
C905637	C905637
C905272	C905272
C901310	C902471

Trunk Positioning

- | | | |
|--|------------------------------|------------------------------|
| <input type="checkbox"/> H-Harness with Padded Covers (Requires 3-Point Pelvic Belt) | Price
\$107 | HCPCS
E0960 |
| <input type="checkbox"/> Full Torso Support Vest | Price
\$355 | K0108 |
| <input type="checkbox"/> Reducer Seat Insert (Reduces Seat Width 2" & Seat Depth 1") | Price
\$417 | K0108 |

SFT16	SFT18
C903988	C903988
C901308	C902469
C904110	C904110

Pelvic Positioning

- | | | |
|--|------------------------------|------------------------------|
| <input type="checkbox"/> 3-Point Pelvic Belt | Price
\$53 | HCPCS
E0978 |
| <input type="checkbox"/> Position Cushion, Medial Thigh Support/Anti Thrust ^{1*} | Price
\$508 | E2607 |
| <input type="checkbox"/> Align Cushion, Lateral Pelvis/Lateral Thigh Support ^{1*} | Price
\$508 | E2607 |
| <input type="checkbox"/> Incontinence Cushion Liner ¹ (Cushion Required) | Price
\$79 | E2619 |

SFT16	SFT18
C904017	C904017
C903590	C903591
C903594	C903595
C902324	C902324

¹Manufacturer Installation Required.

*Solid Seat Insert included with Cushion.

Foot & Leg Positioning

	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Medial Thigh Support, Abductor	\$217	E0957
<input type="checkbox"/> Lateral Thigh Support, Adductor	\$217	E0956
<input type="checkbox"/> Elevating Legrests with Swing-Away Footplates	\$363	E0990
<input type="checkbox"/> 3" Legrest Extension*^	\$74	
<input type="checkbox"/> Calf Panel	\$67	E0995
<input type="checkbox"/> Angle Adjustable Footplates^	\$240	K0040
Foot Positioners - Select Size Below:		
<input type="checkbox"/> Small 7" - 10"	\$170	K0108
<input type="checkbox"/> Medium 8" - 11"	\$170	K0108
<input type="checkbox"/> Large 9" - 12"	\$170	K0108
<input type="checkbox"/> Footplate Securement Strap	\$34	

<u>SFT16</u>	<u>SFT18</u>
C901297	C902458
C901298	C902459
C901299	C902460
C901324	C902478
C904580	C904580
C903271	C903271
C903415	C903415
C903416	C903416
C903417	C903417
C903584	C903584

***3" Legrest Extension interferes with caster swivel when seat is in upright position.
 Seat must be tilted back at least 5 degrees OR Legrests must be elevated to allow caster swivel.
 ^Elevating Legrests with Swing-Away Footplates required.**

Tire Options

	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> 12.5" x 2" Rear Solid Knobby Tires	\$217	
<input type="checkbox"/> 7.5" x 2" Front, 12.5" x 2" Rear Pneumatic Knobby Tires	\$325	
<input type="checkbox"/> 12.5" x 2" Rear Solid Knobby Tires	\$235	

<u>SFT16</u>	<u>SFT18</u>
C904117	
C904118	
	C903420

Accessories

	<u>Price</u>	<u>HCPCS</u>
<input type="checkbox"/> Headrest Cover (Canopy)	\$267	K0108
<input type="checkbox"/> Extended Headrest Cover (Canopy)	\$292	K0108
<input type="checkbox"/> Upper Extremity Support Surface Tray (Clear Tray)^	\$131	E0950
<input type="checkbox"/> Adj. Removeable Hardware (Required for Clear Tray)	\$266	E1028
<input type="checkbox"/> Solid Side Panels - (N/A w/ Recline Option)	\$284	K0108
<input type="checkbox"/> Wheelchair Travel Bag	\$304	
<input type="checkbox"/> Utility Bag	\$46	
<input type="checkbox"/> Mosquito Net (Canopy Required)	\$164	

<u>SFT16</u>	<u>SFT18</u>
C901313	C902475
C904734	C904735
C903403	C903404
C904519	C904519
C904392	C904393
C902022	C902022
C904531	C904531
C904656	C904656

^Height Adjustable Flip-Up Armrest Required

The Safari transit option includes four manufacturer installed wheelchair transportation anchors. The Safari transit model has been crash tested and conforms to RESNA WC4:2017 section 19 and ISO 7176-19:2022 standards.

HCPCS codes are provided for your convenience and are NOT to be considered billing or legal advice. Providers are responsible for insuring that appropriate billing codes are used when submitting claims. Proper use of HCPCS codes does not insure coverage or payment. Contact the appropriate payer for coverage information.

Standard Domestic Shipping Rate: \$40 • For expedited shipping rates, please call 1-844-876-6245 or email convoidsales.us@etac.com

Prices in USD & subject to change, ex works Torrance, California.

Effective: 1/20/2025

Convoid
 2830 California Street
 Torrance, CA 90503
 1 844-US-MOBILITY (844-876-6245)
www.etac.us.com

OF0600-ENG-REV.00

