



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

# Mustang

Company:		Ship To:	
Account Number:	Date:	Address:	
Quote:	PO Number:	City:	Province:
Clinician / Sales Rep:		Postal Code:	Phone:
Marked For:		Email/Fax:	
Associated Quote/Order#:			Effective: 1/20/2025

## Mustang

mm (inch), kg (lb)


	 1	 2	 3	 4
 Overall width 4" or 8" wheels	520 mm (20½")	580 mm (22¾")	650 mm (25½")	700 mm (27½")
 Inside width	380 mm (15")	450 mm (17¾")	500 mm (19½")	560 mm (22")
 Length with 4" wheels	750 mm (29½")	830 mm (32¾")	-	-
 Length with 8" wheels	820 mm (32¼")	900 mm (35½")	1000 mm (39¼")	1100 mm (43¼")
 Height with 4" wheels	400-650 mm (15¾-25½")	600-850 mm (23½-33½")	-	-
 Height with 8" wheels	510-650 mm (20-25½")	600-850 mm (23½-33½")	700-1050 mm (27½-41¼")	1030-1250 mm (40½-49¼")
 Seat height with 4" wheels PU seat (accessory)	170-350 mm (6¾-13¾")	240-410 mm (9½-16¼")	-	-
 Seat height with 8" wheels PU seat (accessory)	260-440 mm (10¼-17¼")	330-500 mm (13-19¾")	450-680 mm (17¾-26¾")	620-840 mm (24½-33")
 Seat height with 4" wheels Comfort seat (accessory)	170-430 mm (6¾-17")	240-490 mm (9½-19¼")	-	-
 Seat height with 8" wheels Comfort seat (accessory)	260-520 mm (10¼-20½")	330-580 mm (13-22¾")	450-760 mm (17¾-30")	620-920 mm (24½-36¼")
 Width between handles (handles = accessory)	165-320 mm (6¼-12½")	165-320 mm (6¼-12½")	140-400 mm (6½-15¾")	140-400 mm (6½-15¾")
 Angle, center bar (a1)	30°	30°	25°	25°
 Minimum turning radius	360 mm (14¼")	550 mm (21¾")	600 mm (23½")	650 mm (26½")
 User height, floor to axil	400-650 mm (15¾-25½")	600-850 mm (23½-33½")	700-1050 mm (27½-41¼")	870-1300 mm (34¼-51¼")



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

## Mustang Frame Options

### Mustang, Standard Package

	<b>Mustang, Standard Package</b>	
	Mustang, Yellow, 4" Casters, Size 1 (869011)	\$2633
	Mustang, Yellow, 8" Casters, Size 1 (869015)	\$2633
	Mustang, Yellow, 8" Casters Slim Line Wheels, Size 1 (869015-SL)	\$2633
	Mustang, White, 4" Casters, Size 2 (869025-4)	\$2742
	Mustang, White, 8" Casters, Size 2 (869024)	\$2742
	Mustang, White, 8" Casters Slim Line Wheels, Size 2 (869024-SL)	\$2742
	Mustang, Red, 8" Casters, Size 3 (869033)	\$3879
	Mustang, Red, 8" Casters Slim Line Wheels, Size 3 (869033-SL)	\$3879
	Mustang, Anthracite (Grey), 8" Casters, Size 4 (869042)	\$4436
Mustang, Anthracite (Grey), 8" Casters Slim Line Wheels, Size 4 (869042-SL)	\$4436	

## Surrounding Supports and Covers

### Surrounding Chest Support and Covers

*Required Selection*



#### Surrounding Chest Support

Chest Support, 6"- 8.5", Small (869810)	\$663
Chest Support, 9.5" - 12", Medium (869820)	\$730
Chest Support, 11"- 14.75", Large (869830)	\$795

#### Covers

Cover for Chest Support, Small (869810-160)	\$158
Cover for Chest Support, Medium (869820-260)	\$191
Cover for Chest Support, Large (869830-360)	\$204

#### Mounting For Chest Support

Low Mounting Kit for Surrounding Chest Support (869800-KIT)	\$332
---	-------

## Seat Options

### Seats and Slings



#### Sling Seat

Sling Seat for Small Chest Support, W: 7" x L: 9.5", Size 1 (869570)	\$239
Sling Seat for Medium Chest Support, W: 7" x L: 12.75", Size 2 (869571)	\$239



#### Mustang Seat

Mustang Seat, Small, Size 1 and 2 (869110)	\$348
Mustang Seat, Large, Size 3 and 4 (869120)	\$368
Mustang Comfort Seat, Small, Size 1 and 2 (869115)	\$504
Mustang Comfort Seat, Large, Size 3 and 4 (869125)	\$525

## Additional Support Options

### Supports



#### Pelvic Supports

*Mounts on center rod*

Surrounding Pelvic Supports, Small (869610)	\$647
Surrounding Pelvic Supports, Medium (869620)	\$667
Surrounding Pelvic Supports, Large (869630)	\$688

#### Extra Pad

Extra Pad for Chest/Surrounding Pelvic Supports (869225)	\$83
--	------

#### Sacrum Support

*Requires Mustang Seat*

Sacrum Support, Small (869250)	\$219
Sacrum Support, Large (869255)	\$299

#### Hip Supports

*Requires Mustang Seat*

Hip Support on Straight Bar, Size 1 and 2 (869221)	\$525
Hip Support on Straight Bar, Size 3 and 4 (869231)	\$525
Hip Support on Curved Bar, Size 1 and 2 (869222)	\$525
Hip Support on Curved Bar, Size 3 and 4 (869232)	\$525



#### Abdominal Support

Abdominal Support (869260)	\$260
----------------------------	-------



# Mustang

Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

## Armrest Options

### Armrests

#### Armrests

Armrest for Anterior Use W/ Handles, Small (869350)	\$558
Armrest for Anterior Use W/ Handles, Large (869355)	\$688

#### Arm Prompts

Arm Prompts for Anterior Use, Small (869310-1)	\$666
Arm Prompts for Anterior Use, Medium (869320-1)	\$666
Arm Prompts for Anterior Use, Large (869325-1)	\$666
Arm Prompts for Posterior Use, Small (869330-1)	\$692
Arm Prompts for Posterior Use, Medium (869340-1)	\$692
Arm Prompts for Posterior Use, Large (869345-1)	\$692



#### Padded Forearm Positioner (Each)

*If wanting a pair of straps, quantity two needs to be selected*

Padded Forearm Positioner Strap (869370)	\$93
--	------

#### Decelerating Handbrake

Decelerating Handbrake, Size 1 (869346-1)	\$1983
Decelerating Handbrake, Size 2 (869346-2)	\$1983
Decelerating Handbrake, Size 3 (869346-2)	\$1983
Decelerating Handbrake, Size 4 (869346-4)	\$1983

## Headrest and Hardware

### Head Support



#### Headrest

Anatomic Head Support, Small (5099001)	\$348
Anatomic Head Support, Medium (5099002)	\$348
Anatomic Head Support, Large (5099003)	\$348

#### Headrest Cover

Anatomic Head Support Cover, Small (9923-83)	N/C
Anatomic Head Support Cover, Medium	N/C
Anatomic Head Support Cover, Large	N/C

#### Headrest Mounting Hardware

Headrest Mounting Hardware, Size 1 and 2, Short (98619-25)	\$314
Headrest Mounting Hardware, Size 3 and 4, Long (98619-26)	\$314

#### Headband System

Headband System (99500)	\$314
-------------------------	-------



#### Connecting Rod

*Required for use of headrest*

Connecting Rod Between Hip and Head Support (869226)	\$77
--	------

## Accessories

### Leg Separator and Ankle Prompts



#### Leg Separators

Leg Separator for Size 1 w/ 4" Wheels (869410)	\$395
Leg Separator for Size 1 w/ 8" Wheels (869410-8)	\$395
Leg Separator for Size 2 w/ 8" Wheels (869420)	\$395
Leg Separator for Size 2 w/ 4" wheels and Size 3 (869430)	\$395
Leg Separator for Size 4 (869440)	\$477



#### Leg Pad

Leg Pad (869470)	\$185
------------------	-------



#### Ankle Prompts

Ankle Prompts for Size 1 (869450)	\$239
Ankle Prompts for Size 2 w/ 8" Wheels (869455)	\$239
Ankle Prompts for Size 2 w/ 4" Wheels (869455-4)	\$239
Ankle Prompts for Size 3 (869460)	\$259
Ankle Prompts for Size 4 (869465)	\$259



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com

## Accessories



### Attendant Guide Bar

Attendant Guide Bar, Size 1 and 2 (869510) \_\_\_\_\_ \$273



### Storage Pouch

Storage Pouch, Size 1 and 2 (869560) \_\_\_\_\_ \$132



### Storage Bag

Storage Bag, 3 and 4 (869565) \_\_\_\_\_ \$265