



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Mustang

Company:		Ship To:	
Account Number:	Date:	Address:	
Quote:	PO Number:	City:	State:
ATP:		Zip:	Phone:
Marked For:		Email/Fax:	
Associated Quote/Order#:			Effective: 1/20/2025

Mustang

mm (inch), kg (lb)

	 1	 2	 3	 4
 Overall width 4" or 8" wheels	520 mm (20½")	580 mm (22¾")	650 mm (25½")	700 mm (27½")
 Inside width	380 mm (15")	450 mm (17¾")	500 mm (19½")	560 mm (22")
 Length with 4" wheels	750 mm (29½")	830 mm (32¾")	-	-
 Length with 8" wheels	820 mm (32¼")	900 mm (35½")	1000 mm (39¼")	1100 mm (43¼")
 Height with 4" wheels	400-650 mm (15¾-25½")	600-850 mm (23½-33½")	-	-
 Height with 8" wheels	510-650 mm (20-25½")	600-850 mm (23½-33½")	700-1050 mm (27½-41¼")	1030-1250 mm (40½-49¼")
 Seat height with 4" wheels PU seat (accessory)	170-350 mm (6¾-13¾")	240-410 mm (9½-16¼")	-	-
 Seat height with 8" wheels PU seat (accessory)	260-440 mm (10¼-17¼")	330-500 mm (13-19¾")	450-680 mm (17¾-26¾")	620-840 mm (24½-33")
 Seat height with 4" wheels Comfort seat (accessory)	170-430 mm (6¾-17")	240-490 mm (9½-19¼")	-	-
 Seat height with 8" wheels Comfort seat (accessory)	260-520 mm (10¼-20½")	330-580 mm (13-22¾")	450-760 mm (17¾-30")	620-920 mm (24½-36¼")
 Width between handles (handles = accessory)	165-320 mm (6¼-12½")	165-320 mm (6¼-12½")	140-400 mm (6½-15¾")	140-400 mm (6½-15¾")
 Angle, center bar (a1)	30°	30°	25°	25°
 Minimum turning radius	360 mm (14¼")	550 mm (21¾")	600 mm (23½")	650 mm (26½")
 User height, floor to axil	400-650 mm (15¾-25½")	600-850 mm (23½-33½")	700-1050 mm (27½-41¼")	870-1300 mm (34¼-51¼")

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.
 Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.
DCN1272.1 Pricing and specifications subject to change without notice.



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Mustang

Mustang Frame Options

Mustang, Standard Package



Mustang, Standard Package

Mustang, Yellow, 4" Casters, Size 1 (869011)	\$2112
Mustang, Yellow, 8" Casters, Size 1 (869015)	\$2112
Mustang, Yellow, 8" Casters Slim Line Wheels, Size 1 (869015-SL)	\$2112
Mustang, White, 4" Casters, Size 2 (869025-4)	\$2199
Mustang, White, 8" Casters, Size 2 (869024)	\$2199
Mustang, White, 8" Casters Slim Line Wheels, Size 2 (869024-SL)	\$2199
Mustang, Red, 8" Casters, Size 3 (869033)	\$3109
Mustang, Red, 8" Casters Slim Line Wheels, Size 3 (869033-SL)	\$3109
Mustang, Anthracite (Grey), 8" Casters, Size 4 (869042)	\$3557
Mustang, Anthracite (Grey), 8" Casters Slim Line Wheels, Size 4 (869042-SL)	\$3557

Surrounding Supports and Covers

Surrounding Chest Support and Covers

Required Selection



Surrounding Chest Support

Chest Support, 6"- 8.5", Small (869810)	\$510
Chest Support, 9.5" - 12", Medium (869820)	\$561
Chest Support, 11"- 14.75", Large (869830)	\$612

Covers

Cover for Chest Support, Small (869810-160)	\$126
Cover for Chest Support, Medium (869820-260)	\$155
Cover for Chest Support, Large (869830-360)	\$165

Mounting For Chest Support

Low Mounting Kit for Surrounding Chest Support (869800-KIT)	\$255
---	-------

Seat Options

Seats and Slings



Sling Seat

Sling Seat for Small Chest Support, W: 7" x L: 9.5", Size 1 (869570)	\$191
Sling Seat for Medium Chest Support, W: 7" x L: 12.75", Size 2 (869571)	\$191



Mustang Seat

Mustang Seat, Small, Size 1 and 2 (869110)	\$279
Mustang Seat, Large, Size 3 and 4 (869120)	\$295
Mustang Comfort Seat, Small, Size 1 and 2 (869115)	\$405
Mustang Comfort Seat, Large, Size 3 and 4 (869125)	\$421

Additional Support Options

Supports



Pelvic Supports

Mounts on center rod

Surrounding Pelvic Supports, Small (869610)	\$520
Surrounding Pelvic Supports, Medium (869620)	\$535
Surrounding Pelvic Supports, Large (869630)	\$552

Extra Pad

Extra Pad for Chest/Surrounding Pelvic Supports (869225)	\$68
--	------

Sacrum Support

Requires Mustang Seat

Sacrum Support, Small (869250)	\$175
Sacrum Support, Large (869255)	\$241

Hip Supports

Requires Mustang Seat

Hip Support on Straight Bar, Size 1 and 2 (869221)	\$421
Hip Support on Straight Bar, Size 3 and 4 (869231)	\$421
Hip Support on Curved Bar, Size 1 and 2 (869222)	\$421
Hip Support on Curved Bar, Size 3 and 4 (869232)	\$421



Abdominal Support

Abdominal Support (869260)	\$200
----------------------------	-------



Phone (800) 981-1540 | Fax (715) 254-0996
 sales@kimobility.com | www.kimobility.com

Mustang

Armrest Options

Armrests

Armrests

Armrest for Anterior Use W/ Handles, Small (869350)	\$447
Armrest for Anterior Use W/ Handles, Large (869355)	\$552

Arm Prompts

Arm Prompts for Anterior Use, Small (869310-1)	\$534
Arm Prompts for Anterior Use, Medium (869320-1)	\$534
Arm Prompts for Anterior Use, Large (869325-1)	\$534
Arm Prompts for Posterior Use, Small (869330-1)	\$555
Arm Prompts for Posterior Use, Medium (869340-1)	\$555
Arm Prompts for Posterior Use, Large (869345-1)	\$555



Padded Forearm Positioner (Each)

If wanting a pair of straps, quantity two needs to be selected

Padded Forearm Positioner Strap (869370)	\$75
--	------

Decelerating Handbrake

Decelerating Handbrake, Size 1 (869346-1)	\$1525
Decelerating Handbrake, Size 2 (869346-2)	\$1525
Decelerating Handbrake, Size 3 (869346-2)	\$1525
Decelerating Handbrake, Size 4 (869346-4)	\$1525

Headrest and Hardware

Head Support



Headrest

Anatomic Head Support, Small (5099001)	\$279
Anatomic Head Support, Medium (5099002)	\$279
Anatomic Head Support, Large (5099003)	\$279

Headrest Cover

Anatomic Head Support Cover, Small (9923-83)	N/C
Anatomic Head Support Cover, Medium	N/C
Anatomic Head Support Cover, Large	N/C

Headrest Mounting Hardware

Headrest Mounting Hardware, Size 1 and 2, Short (98619-25)	\$251
Headrest Mounting Hardware, Size 3 and 4, Long (98619-26)	\$251

Headband System

Headband System (99500)	\$251
-------------------------	-------



Connecting Rod

Required for use of headrest

Connecting Rod Between Hip and Head Support (869226)	\$61
--	------

Accessories

Leg Separator and Ankle Prompts



Leg Separators

Leg Separator for Size 1 w/ 4" Wheels (869410)	\$317
Leg Separator for Size 1 w/ 8" Wheels (869410-8)	\$317
Leg Separator for Size 2 w/ 8" Wheels (869420)	\$317
Leg Separator for Size 2 w/ 4" wheels and Size 3 (869430)	\$317
Leg Separator for Size 4 (869440)	\$382



Leg Pad

Leg Pad (869470)	\$147
------------------	-------



Ankle Prompts

Ankle Prompts for Size 1 (869450)	\$191
Ankle Prompts for Size 2 w/ 8" Wheels (869455)	\$191
Ankle Prompts for Size 2 w/ 4" Wheels (869455-4)	\$191
Ankle Prompts for Size 3 (869460)	\$208
Ankle Prompts for Size 4 (869465)	\$208



Phone (800) 981-1540 | Fax (715) 254-0996
sales@kimobility.com | www.kimobility.com

Mustang

Accessories



Attendant Guide Bar

Attendant Guide Bar, Size 1 and 2 (869510) _____ \$219



Storage Pouch

Storage Pouch, Size 1 and 2 (869560) _____ \$107



Storage Bag

Storage Bag, 3 and 4 (869565) _____ \$214