



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

# Mustang

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 9/5/2024

## Mustang

mm (inch), kg (lb)

|   |  1 |  2 |  3 |  4 |
|---|---|---|---|---|
|  Overall width 4" or 8" wheels   | 520 mm<br>(20½")  | 580 mm<br>(22¾")  | 650 mm<br>(25½")  | 700 mm<br>(27½")  |
|  Inside width  | 380 mm<br>(15")   | 450 mm<br>(17¾")  | 500 mm<br>(19½")  | 560 mm<br>(22")   |
|  Length with 4" wheels                                    | 750 mm<br>(29½")  | 830 mm<br>(32¾")  | -   | -   |
|  Length with 8" wheels                                  | 820 mm<br>(32¼")  | 900 mm<br>(35½")  | 1000 mm<br>(39¼")   | 1100 mm<br>(43¼")   |
|  Height with 4" wheels                                  | 400-650 mm<br>(15¾-25½")  | 600-850 mm<br>(23½-33½")  | -   | -   |
|  Height with 8" wheels                                  | 510-650 mm<br>(20-25½")   | 600-850 mm<br>(23½-33½")  | 700-1050 mm<br>(27½-41¼")   | 1030-1250 mm<br>(40½-49¼")  |
|  Seat height with 4" wheels<br>PU seat (accessory)      | 170-350 mm<br>(6¾-13¾")   | 240-410 mm<br>(9½-16¼")   | -   | -   |
|  Seat height with 8" wheels<br>PU seat (accessory)      | 260-440 mm<br>(10¼-17¼")  | 330-500 mm<br>(13-19¾")   | 450-680 mm<br>(17¾-26¾")  | 620-840 mm<br>(24½-33")   |
|  Seat height with 4" wheels<br>Comfort seat (accessory) | 170-430 mm<br>(6¾-17")  | 240-490 mm<br>(9½-19¼")   | -   | -   |
|  Seat height with 8" wheels<br>Comfort seat (accessory) | 260-520 mm<br>(10¼-20½")  | 330-580 mm<br>(13-22¾")   | 450-760 mm<br>(17¾-30")   | 620-920 mm<br>(24½-36¼")  |
|  Width between handles (handles = accessory)                               | 165-320 mm<br>(6¼-12½")   | 165-320 mm<br>(6¼-12½")   | 140-400 mm<br>(6½-15¾")   | 140-400 mm<br>(6½-15¾")   |
|  Angle, center bar (a1)  | 30°   | 30°   | 25°   | 25°   |
|  Minimum turning radius  | 360 mm<br>(14¼")  | 550 mm<br>(21¾")  | 600 mm<br>(23½")  | 650 mm<br>(26½")  |
|  User height, floor to axil  | 400-650 mm<br>(15¾-25½")  | 600-850 mm<br>(23½-33½")  | 700-1050 mm<br>(27½-41¼")   | 870-1300 mm<br>(34¼-51¼")   |

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.


DCN1272.0 Pricing and specifications subject to change without notice.



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

# Mustang

## Mustang, Standard Package



### Mustang, Standard Package

|   |        |
|---|--------|
| Mustang, Yellow, 4" Casters, Size 1 (869011)                                | \$2070 |
| Mustang, Yellow, 8" Casters, Size 1 (869015)                                | \$2070 |
| Mustang, Yellow, 8" Casters Slim Line Wheels, Size 1 (869015-SL)            | \$2070 |
| Mustang, White, 4" Casters, Size 2 (869025-4)                               | \$2155 |
| Mustang, White, 8" Casters, Size 2 (869024)                                 | \$2155 |
| Mustang, White, 8" Casters Slim Line Wheels, Size 2 (869024-SL)             | \$2155 |
| Mustang, Red, 8" Casters, Size 3 (869033)                                   | \$3048 |
| Mustang, Red, 8" Casters Slim Line Wheels, Size 3 (869033-SL)               | \$3048 |
| Mustang, Anthracite (Grey), 8" Casters, Size 4 (869042)                     | \$3487 |
| Mustang, Anthracite (Grey), 8" Casters Slim Line Wheels, Size 4 (869042-SL) | \$3487 |

## Surrounding Chest Support and Covers



### Surrounding Chest Support

|   |       |
|---|-------|
| Chest Support, 6" - 8.5", Small (869810)    | \$500 |
| Chest Support, 9.5" - 12", Medium (869820)  | \$550 |
| Chest Support, 11" - 14.75", Large (869830) | \$600 |




### Cover for Surrounding Chest Support

|  |       |
|--|-------|
| Cover for Chest Support, Small (869810-160)  | \$123 |
| Cover for Chest Support, Medium (869820-260) | \$151 |
| Cover for Chest Support, Large (869830-360)  | \$161 |

### Mounting For Chest Support


|   |       |
|---|-------|
| Low Mounting Kit for Surrounding Chest Support (869800-KIT) | \$250 |
|---|-------|

## Seats and Slings



### Sling Seat

|   |       |
|---|-------|
| Sling Seat for Small Chest Support, W: 7" x L: 9.5", Size 1 (869570)    | \$187 |
| Sling Seat for Medium Chest Support, W: 7" x L: 12.75", Size 2 (869571) | \$187 |



### Mustang Seat


|  |       |
|--|-------|
| Mustang Seat, Small, Size 1 and 2 (869110)         | \$273 |
| Mustang Seat, Large, Size 3 and 4 (869120)         | \$289 |
| Mustang Comfort Seat, Small, Size 1 and 2 (869115) | \$397 |
| Mustang Comfort Seat, Large, Size 3 and 4 (869125) | \$412 |



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com


# Mustang

## Supports



### Surrounding Pelvic Supports


Surrounding Pelvic Supports, Small (869610) \$509  
 Surrounding Pelvic Supports, Medium (869620) \$524  
 Surrounding Pelvic Supports, Large (869630) \$541



### Hip Support on Straight Bar


*Requires mustang seat or comfort seat*

Hip Support on Straight Bar, Size 1 and 2 (869221) \$412  
 Hip Support on Straight Bar, Size 3 and 4 (869231) \$412



### Extra Pad for Chest/Surrounding Pelvic Supports


Extra Pad for Chest/Surrounding Pelvic Supports (869225) \$66



### Sacrum Support

*Requires mustang seat or comfort seat*


Sacrum Support, Small (869250) \$171  
 Sacrum Support, Large (869255) \$236



### Hip Support on Curved Bar

*Requires mustang seat or comfort seat*

Hip Support on Curved Bar, Size 1 and 2 (869222) \$412  
 Hip Support on Curved Bar, Size 3 and 4 (869232) \$412



### Abdominal Support


Abdominal Support (869260) \$196

## Armrests



### Armrest for Anterior Use W/ Handles

Armrest for Anterior Use W/ Handles, Small (869350) \$438  
 Armrest for Anterior Use W/ Handles, Large (869355) \$541



### Padded Forearm Positioner (Each)

*Sold as an each. Please select two if needing for two arms.*

Padded Forearm Positioner Strap (869370) \$73



### Arm Prompts for Anterior Use

Arm Prompts for Anterior Use, Small (869310-1) \$523  
 Arm Prompts for Anterior Use, Medium (869320-1) \$523  
 Arm Prompts for Anterior Use, Large (869325-1) \$523


### Decelerating Handbrake

Decelerating Handbrake, Size 1 (869346-1) \$1495  
 Decelerating Handbrake, Size 2 (869346-2) \$1495  
 Decelerating Handbrake, Size 3 (869346-2) \$1495  
 Decelerating Handbrake, Size 4 (869346-4) \$1495

### Arm Prompts for Posterior Use


Arm Prompts for Posterior Use, Small (869330-1) \$544  
 Arm Prompts for Posterior Use, Medium (869340-1) \$544  
 Arm Prompts for Posterior Use, Large (869345-1) \$544

## Head Support



### Headrest Mounting Hardware

Headrest Mounting Hardware, Size 1 and 2, Short (98619-25) \$246  
 Headrest Mounting Hardware, Size 3 and 4, Long (98619-26) \$246



### Headrest

Anatomic Head Support, Small (5099001) \$273  
 Anatomic Head Support, Medium (5099002) \$273  
 Anatomic Head Support, Large (5099003) \$273

### Headrest Cover

*Required when picking an Anatomic headrest*

Anatomic Head Support Cover, Small (9923-83) N/C  
 Anatomic Head Support Cover, Medium N/C  
 Anatomic Head Support Cover, Large N/C

### Headband System

Headband System (99500) \$246



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com

# Mustang

## Leg Separator and Ankle Prompts

### Leg Separators



- Leg Separator for Size 1 w/ 4" Wheels (869410) \_\_\_\_\_ \$310
- Leg Separator for Size 1 w/ 8" Wheels (869410-8) \_\_\_\_\_ \$310
- Leg Separator for Size 2 and 3 w/ 4" Wheels (869430) \_\_\_\_\_ \$310
- Leg Separator for Size 2 w/ 8" Wheels (869420) \_\_\_\_\_ \$310
- Leg Separator for Size 4 (869440) \_\_\_\_\_ \$374

### Ankle Prompts



- Ankle Prompts for Size 1 (869450) \_\_\_\_\_ \$187
- Ankle Prompts for Size 2 w/ 8" Wheels (869455) \_\_\_\_\_ \$187
- Ankle Prompts for Size 2 w/ 4" Wheels (869455-4) \_\_\_\_\_ \$187
- Ankle Prompts for Size 3 (869460) \_\_\_\_\_ \$203
- Ankle Prompts for Size 4 (869465) \_\_\_\_\_ \$203

### Leg Pad



- Leg Pad (869470) \_\_\_\_\_ \$144

## Accessories

### Attendant Guide Bar

- Attendant Guide Bar, Size 1 and 2 (869510) \_\_\_\_\_ \$214



### Storage Pouch

- Storage Pouch, Size 1 and 2 (869560) \_\_\_\_\_ \$104

### Storage Bag



- Storage Bag, 3 and 4 (869565) \_\_\_\_\_ \$209