

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

Province:

Clinician / Sales Rep:

Postal Code:

Phone:



Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/20/2025


## Rabbit Up



A (Total Height) \_\_\_\_\_  
B (Armpit Height) \_\_\_\_\_  
C (Top of Pelvic Height) \_\_\_\_\_  
D (Buttock Height) \_\_\_\_\_  
E (Chest Width) \_\_\_\_\_  
F (Pelvic Width) \_\_\_\_\_  
Age \_\_\_\_\_  
Weight \_\_\_\_\_

### Features, incl. in standard price


**Stander Frame**  
Centre Bar  
Brake on Castors  
Angling Spindle




**Mobile Stander Frame**  
Centre Bar  
Driving Wheels  
Angling Spindle  
Closed Pushrim  
0° Camber std. Angle  
Spoke protectors  
Central brake  
Brake on Casters




### Measurements




Size 1



Size 2



Size 3



Size 4

### Product Specs as a standing frame

	Size 1	Size 2	Size 3	Size 4
Overall width (G)	20¼"	21¾"	23¾"	23¾"
Length (H)	27¾"	31¼"	31¼"	36½"
Total mass	19.8 lbs.	23.1 lbs.	24.2 lbs.	29.3 lbs.

### Product Specs as a mobile stander

	w/26" wheels	w/32" wheels	w/32" wheels	w/36" wheels
Overall width (G)	25"	26¾"	28¾"	28¼"
Length (H)	24½"	29½"	29½"	35½"
Total mass	29.75 lbs.	34.1 lbs.	35.25 lbs.	40.75 lbs.



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



# Rabbit Up

	Size 1	Size 2	Size 3	Size 4
<b>Overall measurements</b>				
Width, inside (G1)	12¾"	14¾"	16¾"	16¾"
Overall height (I) *	26 ¼" – 39 ¼"	33 ¾" – 46 ¾"	39" – 51 ½"	47 ¼" – 60 ¼"
Height from foot plate (I1) *	21 ¾" – 29 ¼"	29 ½" – 42 ½"	32 ½" – 46 ¾"	37 ¾" – 56"
Center bar angle (a1)	0 to +30°	0 to +30°	0 to +30°	0 to +30°
Foot plate angle (a2)	+15° to -10°	+15° to -10°	+15° to -10°	+15° to -10°
Camber	8°	6°	6°	4°
Center bar length	19¾"	27½"	31½"	39¼"

<b>User Measurement Guide for Stander</b>				
Recommended user height	30" - 39"	39" - 54"	43" - 58"	48" - 71"
Max armpit height	29"	42"	46"	56"
User weight max	88.1 lbs.	110.2 lbs.	132.2 lbs.	154.3 lbs.
Recommended user age	1 - 3 yr./old	4 - 10 yr./old	5 - 12 yr./old	7 - 18 yr./old

\* Measured in upright position and with the foot plate in the lowest position. For smaller children it is possible to raise the foot plate.

## Rabbit-Up Prone Stander



### Rabbit-up, Stander

- Rabbit-Up, Stander, Yellow, Size 1 (8621220-00) \$1868
- Rabbit-Up, Stander, Blue, Size 2 (8622220-00) \$1915
- Rabbit-Up, Stander, Grey, Size 3 (8623220-00) \$1962
- Rabbit-Up, Stander, Silver, Size 4 (8624220-00) \$2106

## Rabbit, Mobile Prone Stander



- 22" Wheels, Yellow, Size 1 (8621222-21) \$3582
- 24" Wheels, Yellow, Size 1 (8621223-21) \$3582
- 26" Wheels, Yellow, Size 1 (8621224-21) \$3582
- 29" Wheels, Blue, Size 2 (8622225-21) \$4486
- 32" Wheels, Blue, Size 2 (8622226-21) \$4486
- 29" Wheels, Grey, Size 3 (8623225-21) \$4534
- 32" Wheels, Grey, Size 3 (8623226-21) \$4534
- 32" Wheels, Silver, Size 4 (8624226-21) \$4677
- 36" Wheels, Silver, Size 4 (8624227-21) \$4819

## Chest Supports and Pad



### Chest Support

- Chest Support, 11" x 15.5" x 7.5", Size 1 (861704) \$280
- Chest Support, 11" x 20-27" x 7.5", Size 2 and 3 (861702) \$280
- Chest Support, 16" x 25-27", 7.5", Size 4 (861703) \$308




### Surrounding Chest Support

- Surrounding Chest Support, 15.5" - 25.25" Small (861810) \$667
- Surrounding Chest Support, 26.25" - 30.75", Medium (861820) \$694
- Surrounding Chest Support, 34.25" - 44.75" Large (861830) \$735



### Extra Pad for Surrounding Support

- Extra Pad for Surrounding Support, H: 4.75" x W: 2.25" (869225) \$83



### Filler Pad

- Filler Pad, Small (861421) \$113
- Filler Pad, Large (861422) \$113



### Extension for Center Bar

- Extension for Center Bar (861480) \$113

### Cover for Chest Support

- Chest Support, Small (869810-160) \$164
- Chest Support, Medium (869820-260) \$200
- Chest Support, Large (869830-360) \$215

## Trays

### Trays



Tray, Plastic, Grey (862515)	\$525
Tray With Storage Box (862501)	\$538
Plastic Tray with bowl and cover (862525)	\$765

### Tray Accessories

Height Adjustment of Tray (862510)	\$308
<i>Required with Tray with Storage Box</i>	

### Mounting Hardware

*Required Selection for Tray with Bowl*

Removeable Hardware Kit for Low Mount of Tray on Center Bar (862525-Kit)	\$378
--	-------

## Abdominal Support



### Abdominal Supports

Abdominal Support (Ea.) (861411)	\$219
----------------------------------	-------

## Head Support



### Connecting System for Head Support


*Requires Hip Support with pommel*

Connecting System for Head Support, Size 1 and 2 (861424-3)	\$352
Connecting System for Head Support, Size 3 and 4 (861424-4)	\$352

### Headrest Cover


*Required Headrest selection*

Anatomic Head Support Cover, Small (9923-83)	N/C
Anatomic Head Support Cover, Medium (9945-83)	N/C
Anatomic Head Support Cover, Large (9967-83)	N/C



### Headrest


Anatomic Head Support, Small (5099001)	\$348
Anatomic Head Support, Medium (5099002)	\$348
Anatomic Head Support, Large (5099003)	\$348



### Headrest Mounting Hardware


Headrest Mounting Hardware, Short (98619-25)	\$314
Headrest Mounting Hardware, Long (98619-26)	\$314

## Pelvic and Hip Control



### Hip Support w/ Pommel

Hip Support w/ Pommel, Size 1 (861501)	\$688
Hip Support w/ Pommel, Size 2 and 3 (861502)	\$688
Hip Support w/ Pommel, Size 4 (861503)	\$763



### Hip Support, With Rotation


*Front mount to center bar*

Hip Support, With Rotation, Size 2 and 3 (861612)	\$484
Hip Support, With Rotation, Size 4 (861613)	\$525



### Hip Support, Front Mount


Hip Support, Front Mount, Size 2 and 3 (861602)	\$257
Hip Support, Front Mount, Size 4 (861603)	\$362



### Hip Support, Surrounding

*Front mount to center bar*

Hip Support, Surrounding, Size 1 (861610)	\$754
Hip Support, Surrounding, Size 2 and 3 (861620)	\$754
Hip Support, Surrounding, Size 4 (861630)	\$770



### Pommel With Back Support

Pommel With Back Support	\$618
--------------------------	-------

## Knee Supports



### Knee Supports (Pair)

Knee Supports, Size 1 (Pair) (861301)	\$381
Knee Supports, Size 2 and 3 (Pair) (861302)	\$484
Knee Supports, Size 4 (Pair) (861303)	\$484




### Hinged Knee Supports (Pair)

Hinged Knee Supports, Size 1 (Pair) (861311)	\$593
Hinged Knee Supports, Size 2 and 3 (Pair) (861312)	\$593
Hinged Knee Supports, Size 4 (Pair) (861313)	\$593

## Footplates

*Flip-up Footplate required for All Sandals, except the Individual Footrest w/ Hinged Sandals*



### Flip-Up Footplate

Flip-Up Footplate, Size 1 and 2 (861211)	\$299
Flip-Up Footplate, Size 2, 3 and 4 (861212)	\$299



### Sandals

Sandals for Rabbit Up/ Toucan, Small (869490-RT)	\$290
Sandals for Rabbit Up/ Toucan, Medium (869491-RT)	\$304
Sandals for Rabbit Up/ Toucan, Large (869492-RT)	\$314



### Individual Footrest w/ Hinged Sandals (Pair)

Individual Footrest w/ Hinged Sandals, Size 1 (Pair) (861201)	\$538
Individual Footrest w/ Hinged Sandals, Size 2 and 3 (Pair) (861202)	\$538

### Hinged Sandles

Hinged Sandals w/ Straps, Size 1 (Pair) (861221)	\$185
Hinged Sandals w/ Straps, Size 2 and 3 (Pair) (861222)	\$185
Hinged Sandals w/ Straps, Size 4 (Pair) (861223)	\$297




### Heel Stops

Heel Stops, Small (Pair) (861251)	\$77
Heel Stops, Large (Pair) (861252)	\$77

### Formed Footplates (Pair)

Formed Footplates, Size 0 (Pair) (86149)	\$226
Formed Footplates, Size 1 and 2 (Pair) (86147)	\$226
Formed Footplates, Size 3 (Pair) (86148)	\$226



### Footplate Straps

Footplate Straps, Plastic Lock (Pair) (85096)	\$77
Footplate Straps, Velcro (Pair) (85102)	\$70

## Conversion Kits

*Please quote two for a complete set of wheels. Two wheels are required for a full conversion kit.*

### Driving Wheel Conversion Kits (Each)

Driving Wheels 22", Closed Pushrims Size 1 (Each) (134209)	\$788
Driving Wheels 24", Closed Pushrims Size 1 (Each) (134210)	\$788
Driving Wheels 26", Closed Pushrims Size 1, 2 and 3 (Each) (134211)	\$788
Driving Wheels 29", Closed Pushrims Size 2 and 3 (Each) (134212)	\$1214
Driving Wheels 32", Closed Pushrims Size 2, 3 and 4 (Each) (134213)	\$1214
Driving Wheels 36", Closed Pushrims Size 4 (Each) (134214)	\$1214



### Camber Kit Rabbit Up

Camber Kit Rabbit Up 0° (134200)	\$71
Camber Kit Rabbit Up +3/-3° (134201)	\$71
Camber Kit Rabbit Up +6/-6° (134202)	\$71



### Central Brake System

Central Brake System, Size 1 (134215)	\$395
Central Brake System, Size 2 and 3 (134216)	\$395
Central Brake System, Size 4 (134218)	\$395