



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

# Rabbit Up

Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:



Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 1/20/2025


## Rabbit Up



- A (Total Height) \_\_\_\_\_
- B (Armpit Height) \_\_\_\_\_
- C (Top of Pelvic Height) \_\_\_\_\_
- D (Buttock Height) \_\_\_\_\_
- E (Chest Width) \_\_\_\_\_
- F (Pelvic Width) \_\_\_\_\_
- Age \_\_\_\_\_
- Weight \_\_\_\_\_

### Features, incl. in standard price


- Stander Frame**
- Centre Bar
  - Brake on Castors
  - Angling Spindle




- Mobile Stander Frame**
- Centre Bar
  - Driving Wheels
  - Angling Spindle
  - Closed Pushrim
  - 0° Camber std. Angle
  - Spoke protectors
  - Central brake
  - Brake on Casters




### Measurements




Size 1



Size 2



Size 3



Size 4

### Product Specs as a standing frame

	Size 1	Size 2	Size 3	Size 4
Overall width (G)	20 1/4"	21 3/4"	23 3/4"	23 3/4"
Length (H)	27 3/4"	31 1/4"	31 1/4"	36 1/2"
Total mass	19.8 lbs.	23.1 lbs.	24.2 lbs.	29.3 lbs.

### Product Specs as a mobile stander

	w/26" wheels	w/32" wheels	w/32" wheels	w/36" wheels
Overall width (G)	25"	26 3/4"	28 3/4"	28 1/4"
Length (H)	24 1/2"	29 1/2"	29 1/2"	35 1/2"
Total mass	29.75 lbs.	34.1 lbs.	35.25 lbs.	40.75 lbs.

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1257.2 Pricing and specifications subject to change without notice.



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

# Rabbit Up

	Size 1	Size 2	Size 3	Size 4
<b>Overall measurements</b>				
Width, inside (G1)	12¾"	14¾"	16¾"	16¾"
Overall height (I) *	26 ¼" – 39 ¼"	33 ¾" – 46 ¾"	39" – 51 ½"	47 ¼" – 60 ¼"
Height from foot plate (I1) *	21 ¾" – 29 ¼"	29 ½" – 42 ½"	32 ½" – 46 ¾"	37 ¾" – 56"
Center bar angle (a1)	0 to +30°	0 to +30°	0 to +30°	0 to +30°
Foot plate angle (a2)	+15° to -10°	+15° to -10°	+15° to -10°	+15° to -10°
Camber	8°	6°	6°	4°
Center bar length	19¾"	27½"	31½"	39¼"

### User Measurement Guide for Stander

Recommended user height	30" - 39"	39" - 54"	43" - 58"	48" - 71"
Max armpit height	29"	42"	46"	56"
User weight max	88.1 lbs.	110.2 lbs.	132.2 lbs.	154.3 lbs.
Recommended user age	1 - 3 yr./old	4 - 10 yr./old	5 - 12 yr./old	7 - 18 yr./old

\* Measured in upright position and with the foot plate in the lowest position. For smaller children it is possible to raise the foot plate.

## Rabbit-Up Prone Stander



### Rabbit-up, Stander

- Rabbit-Up, Stander, Yellow, Size 1 (8621220-00) \$1498
- Rabbit-Up, Stander, Blue, Size 2 (8622220-00) \$1536
- Rabbit-Up, Stander, Grey, Size 3 (8623220-00) \$1574
- Rabbit-Up, Stander, Silver, Size 4 (8624220-00) \$1689

## Rabbit, Mobile Prone Stander



- 22" Wheels, Yellow, Size 1 (8621222-21) \$2872
- 24" Wheels, Yellow, Size 1 (8621223-21) \$2872
- 26" Wheels, Yellow, Size 1 (8621224-21) \$2872
- 29" Wheels, Blue, Size 2 (8622225-21) \$3597
- 32" Wheels, Blue, Size 2 (8622226-21) \$3597
- 29" Wheels, Grey, Size 3 (8623225-21) \$3636
- 32" Wheels, Grey, Size 3 (8623226-21) \$3636
- 32" Wheels, Silver, Size 4 (8624226-21) \$3751
- 36" Wheels, Silver, Size 4 (8624227-21) \$3866

## Chest Supports and Pad



### Chest Support

- Chest Support, 11" x 15.5" x 7.5", Size 1 (861704) \$224
- Chest Support, 11" x 20-27" x 7.5", Size 2 and 3 (861702) \$224
- Chest Support, 16" x 25-27", 7.5", Size 4 (861703) \$246



### Surrounding Chest Support

- Surrounding Chest Support, 15.5" - 25.25" Small (861810) \$535
- Surrounding Chest Support, 26.25" - 30.75", Medium (861820) \$557
- Surrounding Chest Support, 34.25" - 44.75" Large (861830) \$590



### Extra Pad for Surrounding Support

- Extra Pad for Surrounding Support, H: 4.75" x W: 2.25" (869225) \$68



### Filler Pad

- Filler Pad, Small (861421) \$89
- Filler Pad, Large (861422) \$89



### Extension for Center Bar

- Extension for Center Bar (861480) \$89

### Cover for Chest Support

- Chest Support, Small (869810-160) \$126
- Chest Support, Medium (869820-260) \$155
- Chest Support, Large (869830-360) \$165



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

# Rabbit Up

## Trays

### Trays



Tray, Plastic, Grey (862515) \_\_\_\_\_ \$421

Tray With Storage Box (862501) \_\_\_\_\_ \$432

Plastic Tray with bowl and cover (862525) \_\_\_\_\_ \$611

### Tray Accessories

Height Adjustment of Tray (862510) \_\_\_\_\_ \$246  
*Required with Tray with Storage Box*

### Mounting Hardware

*Required Selection for Tray with Bowl*

Removeable Hardware Kit for Low Mount of Tray on Center Bar (862525-Kit) \_\_\_\_\_ \$291

## Abdominal Support



### Abdominal Supports

Abdominal Support (Ea.) (861411) \_\_\_\_\_ \$175

## Head Support



### Connecting System for Head Support

Connecting System for Head Support, Size 1 and 2 (861424-3) \_\_\_\_\_ \$282

Connecting System for Head Support, Size 3 and 4 (861424-4) \_\_\_\_\_ \$282


### Headrest Cover

*Required Headrest selection*

Anatomic Head Support Cover, Small (9923-83) \_\_\_\_\_ N/C

Anatomic Head Support Cover, Medium (9945-83) \_\_\_\_\_ N/C

Anatomic Head Support Cover, Large (9967-83) \_\_\_\_\_ N/C




### Headrest

Anatomic Head Support, Small (5099001) \_\_\_\_\_ \$279

Anatomic Head Support, Medium (5099002) \_\_\_\_\_ \$279

Anatomic Head Support, Large (5099003) \_\_\_\_\_ \$279




### Headrest Mounting Hardware

Headrest Mounting Hardware, Short (98619-25) \_\_\_\_\_ \$251

Headrest Mounting Hardware, Long (98619-26) \_\_\_\_\_ \$251

## Pelvic and Hip Control




### Hip Support w/ Pommel

*Front mount to center bar*

Hip Support w/ Pommel, Size 1 (861501) \_\_\_\_\_ \$552

Hip Support w/ Pommel, Size 2 and 3 (861502) \_\_\_\_\_ \$552


Hip Support w/ Pommel, Size 4 (861503) \_\_\_\_\_ \$611



### Hip Support, With Rotation

Hip Support, With Rotation, Size 2 and 3 (861612) \_\_\_\_\_ \$387


Hip Support, With Rotation, Size 4 (861613) \_\_\_\_\_ \$421



### Hip Support, Front Mount

Hip Support, Front Mount, Size 2 and 3 (861602) \_\_\_\_\_ \$290

Hip Support, Front Mount, Size 4 (861603) \_\_\_\_\_ \$320




### Hip Support, Surrounding

*Front mount to center bar*

Hip Support, Surrounding, Size 1 (861610) \_\_\_\_\_ \$606

Hip Support, Surrounding, Size 2 and 3 (861620) \_\_\_\_\_ \$606

Hip Support, Surrounding, Size 4 (861630) \_\_\_\_\_ \$617



### Pommel With Back Support

Pommel With Back Support \_\_\_\_\_ \$476



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com

# Rabbit Up

## Knee Supports



### Knee Supports (Pair)

Knee Supports, Size 1 (Pair) (861301)	\$305
Knee Supports, Size 2 and 3 (Pair) (861302)	\$387
Knee Supports, Size 4 (Pair) (861303)	\$387




### Hinged Knee Supports (Pair)

Hinged Knee Supports, Size 1 (Pair) (861311)	\$476
Hinged Knee Supports, Size 2 and 3 (Pair) (861312)	\$476
Hinged Knee Supports, Size 4 (Pair) (861313)	\$476

## Footplates

*Flip-up Footplate required for All Sandals, except the Individual Footrest w/ Hinged Sandals*



### Flip-Up Footplate

Flip-Up Footplate, Size 1 and 2 (861211)	\$241
Flip-Up Footplate, Size 2, 3 and 4 (861212)	\$241



### Sandals

Sandals for Rabbit Up/ Toucan, Small (869490-RT)	\$233
Sandals for Rabbit Up/ Toucan, Medium (869491-RT)	\$244
Sandals for Rabbit Up/ Toucan, Large (869492-RT)	\$251

### Formed Footplates (Pair)

Formed Footplates, Size 0 (Pair) (86149)	\$180
Formed Footplates, Size 1 and 2 (Pair) (86147)	\$180
Formed Footplates, Size 3 (Pair) (86148)	\$180

### Hinged Sandles

Hinged Sandals w/ Straps, Size 1 (Pair) (861221)	\$147
Hinged Sandals w/ Straps, Size 2 and 3 (Pair) (861222)	\$147
Hinged Sandals w/ Straps, Size 4 (Pair) (861223)	\$257




### Individual Footrest w/ Hinged Sandals (Pair)

Individual Footrest w/ Hinged Sandals, Size 1 (Pair) (861201)	\$432
Individual Footrest w/ Hinged Sandals, Size 2 and 3 (Pair) (861202)	\$432



### Heel Stops

Heel Stops, Small (Pair) (861251)	\$61
Heel Stops, Large (Pair) (861252)	\$61



### Footplate Straps

Footplate Straps, Plastic Lock (Pair) (85096)	\$40
Footplate Straps, Velcro (Pair) (85102)	\$56

## Conversion Kits

*Please quote two for a complete set of wheels. Two wheels are required for a full conversion kit.*

### Driving Wheel Conversion Kits (Each)

Driving Wheels 22", Closed Pushrims Size 1 (Each) (134209)	\$631
Driving Wheels 24", Closed Pushrims Size 1 (Each) (134210)	\$631
Driving Wheels 26", Closed Pushrims Size 1, 2 and 3 (Each) (134211)	\$631
Driving Wheels 29", Closed Pushrims Size 2 and 3 (Each) (134212)	\$974
Driving Wheels 32", Closed Pushrims Size 2, 3 and 4 (Each) (134213)	\$974
Driving Wheels 36", Closed Pushrims Size 4 (Each) (134214)	\$974



### Camber Kit Rabbit Up

Camber Kit Rabbit Up 0° (134200)	\$57
Camber Kit Rabbit Up +3/-3° (134201)	\$57
Camber Kit Rabbit Up +6/-6° (134202)	\$57



### Central Brake System

Central Brake System, Size 1 (134215)	\$317
Central Brake System, Size 2 and 3 (134216)	\$317
Central Brake System, Size 4 (134218)	\$317