




Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



Company:

Ship To:

Account Number:

Date:

Address:

Quote: PO Number:

City:

State:

ATP:

Zip:

Phone:

Marked For:

Email/Fax:

Associated Quote/Order#:

Effective: 2/3/2025

## x:Panda

mm (inch) kg (lb)

	1	2	3	4
Effective seat width	160, 200, 240 mm (6¼", 7¾", 9½")	200, 240, 280 mm (7¾", 9½", 11")	260, 310, 360 mm (10¼", 12¼", 14¼")	340, 400, 460 mm (13½", 15¾", 18")
Width, thighs	70, 90, 110 mm (2¾", 3½", 4¼")	100, 120, 140 mm (4", 4¾", 5½")	120, 145, 170 mm (4¾", 5¾", 6¾")	130, 190, 250 mm (5", 7½", 9¾")
Effective seat depth	160-260 mm (6¼"-10¼")	220-330 mm (8½"-13")	260-390 mm (10¼"-15¼")	390-530 mm (15¼"-20¾")
Backrest height	270-370 mm (10½"-14½")	320-420 mm (12½"-16½")	350-470 mm (13¾"-18½")	500-710 mm (19¾"-28")
Back width, top of backrest	230 mm (9")	280 mm (11")	330 mm (13")	350+50 mm (13¾"+2")
Overall width, both legs in max abduction	310, 350, 390 mm (12¼", 13¾", 15¼")	380, 420, 460 mm (15", 16½", 18")	470, 520, 570 mm (18½", 20½", 22½")	590, 650, 710 mm (23¾", 25½", 28")
Armrest to seat distance	100-165 mm (4" - 6½")	100-165 mm (4" - 6½")	100-165 mm (4" - 6½")	100-165 mm (4" - 6½")
Front location of armrest structure	140 mm (5½")	140 mm (5½")	180 mm (7")	350 mm (13")
Backrest angle	(-5°)-(30°)	(-5°)-(30°)	(-5°)-(30°)	(-5°)-(30°)
Abduction angle, each leg (s1) Size 1 in the smallest position	0°-15° 0°-8°	0°-15°	0°-15°	0°-15°
Maximum user mass	30 kg (66.1 lb)	30 kg (66.1 lb)	60 kg (132.2 lb)	80 kg (176.3 lb)
Maximum user mass at transportation	30 kg (66.1 lb)	30 kg (66.1 lb)	60 kg (132.2 lb)	75 kg (165.3 lb)
Mass device*	6.5 kg (14.3 lb)	8.5 kg (18.7 lb)	11 kg (24.3 lb)	13 kg (28.5 lb)
Max rated load	38.5 kg (84.7 lb)	40.5 kg (89.2 lb)	73 kg (160.9 lb)	95 kg (209.4 lb)
Minimum User height	650 mm (25½")	850 mm (33½")	1000 mm (39½")	1390 mm (54¾")
Maximum User height	1050 mm (41¼")	1280 mm (50½")	1430 mm (56¼")	1800 mm (70¾")

\*seat only (without accessories)

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.

Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.

DCN1296.2 Pricing and specifications subject to change without notice.



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

x:panda

## x:panda Base Options

*Seat and Back must match Base*

### x:panda



#### Size 1 Panda

Size 1 x:panda, Small (8910111-02)	N/C
Size 1 x:panda, Medium (8910112-02)	N/C
Size 1 x:panda, Large (8910113-02)	N/C

#### Size 2 Panda

Size 2 x:panda, Small (8910121-02)	N/C
Size 2 x:panda, Medium (8910122-02)	N/C
Size 2 x:panda, Large (8910123-02)	N/C

#### Size 4 Panda

Size 4 x:panda, Small (8910141-02)	N/C
Size 4 x:panda, Medium (8910142-02)	N/C
Size 4 x:panda, Large (8910143-02)	N/C

#### Size 3 Panda

Size 3 x:panda, Small (8910131-02)	N/C
Size 3 x:panda, Medium (8910132-02)	N/C
Size 3 x:panda, Large (8910133-02)	N/C

## x:Panda Back and Cushions

### X:Panda Back

*Back selection required. Back size must match above base size.*

#### x:panda Back

Size 1 Back x:panda, Small (X8910111-01) (E2291)	\$710
Size 1 Back x:panda, Medium (X8910112-01) (E2291)	\$710
Size 1 Back x:panda, Large (X8910113-01) (E2291)	\$710
Size 2 Back x:panda, Small (X8910121-01) (E2291)	\$710
Size 2 Back x:panda, Medium (X8910122-01) (E2291)	\$710
Size 2 Back x:panda, Large (X8910123-01) (E2291)	\$710
Size 3 Back x:panda, Small (X8910131-01) (E2291)	\$710
Size 3 Back x:panda, Medium (X8910132-01) (E2291)	\$710
Size 3 Back x:panda, Large (X8910133-01) (E2291)	\$710
Size 4 Back x:panda, Small (X8910141-01) (E2291)	\$727
Size 4 Back x:panda, Medium (X8910142-01) (E2291)	\$727
Size 4 Back x:panda, Large (X8910143-01) (E2291)	\$727

### x:Panda Seat

*Seat selection required. Seat size must match above base size.*

#### x:panda Seat

Size 1 Seat x:panda, Small (X8910111-02) (E2292)	\$727
Size 1 Seat x:panda, Medium (X8910112-02) (E2292)	\$727
Size 1 Seat x:panda, Large (X8910113-02) (E2292)	\$727
Size 2 Seat x:panda, Small (X8910121-02) (E2292)	\$727
Size 2 Seat x:panda, Medium (X8910122-02) (E2292)	\$727
Size 2 Seat x:panda, Large (X8910123-02) (E2292)	\$727
Size 3 Seat x:panda, Small (X8910131-02) (E2292)	\$727
Size 3 Seat x:panda, Medium (X8910132-02) (E2292)	\$727
Size 3 Seat x:panda, Large (X8910133-02) (E2292)	\$727
Size 4 Seat x:panda, Small (X8910141-02) (E2292)	\$727
Size 4 Seat x:panda, Medium (X8910142-02) (E2292)	\$727
Size 4 Seat x:panda, Large (X8910143-02) (E2292)	\$727

### Backrest Options

*Must choose a backrest option. Must match base size.*

#### Dynamic Back

Dynamic Back, 100 Newton Gas Spring, Size 1 and 2 (X911100) (K0108)	\$769
Dynamic Back, 150 Newton Gas Spring, Size 3 (X8911150) (K0108)	\$769
Dynamic Back, 350 Newton Gas Spring, Size 4 (X8911350-4)	\$769

#### Non-Dynamic Back

Non-Dynamic Gas Size 1 and 2 Spring (8911001) (E1014)	\$473
Non-Dynamic Gas Spring, Size 3 (8911002) (E1014)	\$473



#### Recline Option

*Required Selection*

Recline, '-5° to 30° (X66-8911) (E1014)	\$473
---	-------

### Seat Cushion Type

*Required selection*


Split Seat Cushion	N/C	Seat Cushion	N/C
--------------------	-----	--------------	-----

*For use with Medial Thigh System or Abduction Supports*

*For use WITHOUT Medial Thigh System or Abduction Supports*



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



### Seat Cushion Upholstery Color

**Required Selection**

Black Cushion Cover	N/C
Blue Cushion Cover	N/C
Red Cushion Cover	N/C
Pink Cushion Cover	N/C
<i>Not available with size 4</i>	

### Back Cushion Upholstery Color

**Required selection**

Blue Upholstery	N/C
Black Upholstery	N/C
Red Upholstery	N/C
Pink Upholstery	N/C
<i>Not available with size 4</i>	

### Replacement Black Wipeable Seat Cover

**Additional Seat cover only**

#### Size 1 Black Wipeable Cover

Seat cover, Black Wipeable, Size 1, Small (8911411-82)	\$125
Seat cover, Black Wipeable, Size 1, Medium (8911412-82)	\$125
Seat cover, Black Wipeable, Size 1, Large (8911413-82)	\$125

#### Size 2 Black Wipeable Cover

Seat cover, Black Wipeable, Size 2, Small (8911421-82)	\$125
Seat cover, Black Wipeable, Size 2, Medium (8911422-82)	\$125
Seat cover, Black Wipeable, Size 2, Large (8911423-82)	\$125

#### Size 3 Black Wipeable Cover

Seat cover, Black Wipeable, Size 3, Small (8911431-82)	\$125
Seat cover, Black Wipeable, Size 3, Medium (8911432-82)	\$125
Seat cover, Black Wipeable, Size 3, Large (8911433-82)	\$125

#### Size 4 Black Wipeable Cover

Seat cover, Black Wipeable, Size 4, Small (8911441-82)	\$295
Seat cover, Black Wipeable, Size 4, Medium (8911441-82)	\$295
Seat cover, Black Wipeable, Size 4, Large (8911443-82)	\$295

### Replacement Black Wipeable Back Cover

**Additional Seat cover only**

#### Size 1 Back Black Wipeable Cover

Back Cover, Black Wipeable, Size 1, Small	\$125
Back Cover, Black Wipeable, Size 1, Medium	\$125
Back Cover, Black Wipeable, Size 1, Large	\$125

#### Size 2 Back Black Wipeable Cover

Back Cover, Black Wipeable, Size 2, Small	\$125
Back Cover, Black Wipeable, Size 2, Medium	\$125
Back Cover, Black Wipeable, Size 2, Large	\$125

#### Size 3 Black Wipeable

Back Cover, Black Wipeable, Size 3, Small	\$125
Back Cover, Black Wipeable, Size 3, Medium	\$125
Back Cover, Black Wipeable, Size 3, Large	\$125

#### Size 4 Black Wipeable Cover

Back Cover, Black Wipeable, Size 4, Small	\$295
Back Cover, Black Wipeable, Size 4, Medium	\$295
Back Cover, Black Wipeable, Size 4, Large	\$295

## Knee and Postural Supports

### Seat Accessories



#### Anti-Thrust Wedge

Anti-Thrust Wedges (Pr.), Size 1 and 2 (8910400-2) (K0108)	\$55
Anti-Thrust Wedges (Pr.), Size 2 and 3 (8910400-4)(K0108)	\$55
Anti-Thrust Wedges (Pr.), Size 4 (8910400-12) (K0108)	\$75

### Medial Thigh System



#### Medial Thigh System

Medial Thigh System, 3.5" x 3.5" x 3.5" (E0957,E1028) (8910408-1)	\$402
Medial Thigh System, 4.75" x 3.5" x 3.5" (E0957,E1028) (8910408-2)	\$402
Medial Thigh System, 6" x 3.5" x 3.5" (E0957,E1028) (8910408-3)	\$402

### Abduction Supports Left and Right

#### Abduction Support w/ Cover

Abduction Supports Left and Right w/ Covers, Size 1 and 2 (8910402-10)	\$237
Abduction Supports Left and Right w/ Covers, Size 3 (8910402-20)	\$237
Abduction Supports Left and Right w/ Covers, Size 4 (8910402-40)	\$237

### Hip Supports Left and Right w/ Covers



#### Hip Supports (Pair)


Hip Supports Left and Right w/ Covers, Size 1 and 2 (E0956) (8910404-10)	\$237
Hip Supports Left and Right w/ Covers, Size 3 (E0956) (8910404-20)	\$237
Hip Supports Left and Right w/ Covers, Size 4 (E0956) (8910404-40)	\$237

### Hip Supports with PU Pad

Hip Support Left and Right w/ PU Pad, 3.25" x 7.5" (E0956) (8910404-45)	\$644
---	-------



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



## Arm Supports

### Armrest



#### Armrest with Cover

Armrest with Cover, Size 1, (Pair) (E0973) (8910615-2)	\$343
Armrest with Cover, Size 2, (Pair) (E0973) (8910616-2)	\$343
Armrest with Cover, Size 3, (Pair) (E0973) (8910617-2)	\$343
Armrest, Single Post, Size 4 (Pair) (E0973) (8910619)	\$343



#### Armrest Kits

*Used for Armrests 1-3*

Daily Adjustment Kit (8910615-11)	\$33
-----------------------------------	------

### Upper Arm Support



#### Upper Arm Support

Fixed Upper Arm Support with Cover, 5" x 4" (8910407-11)(K0108)	\$169
--	-------

## Lateral Support Cushion and Fittings

### Fixed Laterals

#### Fixed Side Supports w/ Covers (Pair)

Fixed Side Supports Left and Right w/ Covers, Type 1 (8910301-1)	\$237
Fixed Side Supports Left and Right w/ Covers, Type 2 (8910301-2)	\$237
Fixed Side Supports Left and Right w/ Covers, Type 3 (8910301-3)	\$237
Fixed Side Supports Left and Right w/ Cover, Type 4 (8910301-4)	\$237
Fixed Side Supports Left and Right w/ Covers, Type 5 (8910301-5)	\$237
Fixed Side Supports Left and Right w/ Covers, Type 6 (8910301-6)	\$237

### Bars

*Required selection for Swing Away Laterals. Quantity 2 required to mount left and right sides.*

#### Straight Bars

Bar, Straight (81932) (Each)	N/C
------------------------------	-----

#### Offset Bars

Bar, Offset, 1" (81968) (Each)	N/C
Bar, Offset, 1.75" (82000) (Each)	N/C
Bar, Offset, 2.5" (81969) (Each)	N/C

### Swing Away Lateral Support

#### Swing Away Lateral Bracket - Left

Swing-Away Lateral Support Bracket, Left (Each) (83900-L) (E1028)	\$267
---	-------

#### Swing Away Lateral Bracket - Right

Swing-Away Lateral Support Bracket, Right (Each) (83900-R) (E1028)	\$267
--	-------

### Lateral Cushions

*Required selection for laterals. Sold in each, select quantity 2 if wanting both left and right side.*

#### Curved Lateral Cushion

Curved Cushion (Each), 3.5" x 7.5" (6175406) (E0956)	\$125
Curved Cushion (Each), 2.5" x 6.75" (6175405) (E0956)	\$125
Curved Cushion (Each), 2.5" x 6" (6175404) (E0956)	\$125

#### Standard Lateral Cushion

Standard Cushion (Each), 4" x 4" (6175402) (E0956)	\$125
Standard Cushion (Each), 2 2/3" x 4.5" (6175401) (E0956)	\$125

#### Protective Cushion

*For lateral bar*

Protective Cushion for Laterals, Black (89150) (Each)	\$110
---	-------



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com

x:panda

## Positioning Accessories

### Positioning Belts

#### Torso Vest



Torso Vest Small - S-P: 10.75" - 12.5" (85705-S) (E0960)	\$237
Torso Vest Medium - S-P: 12.5" - 16" (85705-M) (E0960)	\$237
Torso Vest Large - S-P: 16" - 18.5" (85705-L) (E0960)	\$237
Torso Vest X-Large - S-P: 18.5" - 20.5" (85705-XL) (E0960)	\$237

#### 4 Point Padded Hip Belt



4-Point Padded Hip Belts 1" X-Small - H-H: 6.75" - 12.5" (85725-XS) (E0960)	\$160
4-Point Padded Hip Belts 1" Small - H-H: 8.25" - 13.75" (85725-S) (E0960)	\$160
4-Point Padded Hip Belts 1" Medium - H-H: 9.75" - 16.25" (85725-M) (E0960)	\$160
4-Point Padded Hip Belts 1" Large - H-H: 12.5" - 17.25" (85725-L) (E0960)	\$160

#### Cross Vest



Cross Vest X-Small - S-P: 7.75" - 12.5" (Each) (85720-XS) (E0960)	\$225
Cross Vest Small - S-P: 10.25" - 13.5" (Each) (85720-S) (E0960)	\$225
Cross Vest Medium - S-P: 13.5" - 17.75" (Each) (85720-M) (E0960)	\$225
Cross Vest Large - S-P: 17.75" - 20.75" (Each) (85720-L) (E0960)	\$225
Cross Vest X-Large - S-P: 20.75" - 24.5" (Each) (85720-XL) (E0960)	\$225

#### 2 Point Pelvic Belt



2-Point Pelvic Belt Black - L: 13-14.5" (85424-1) (E0978)	\$115
---	-------

#### Pelvic Harness



Pelvic Harness, 0 - A:6.5" x B:7.5" (89108-0) (K0108)	\$290
Pelvic Harness, 1 - A:8.5" x B:8.75" (89108-1) (K0108)	\$290
Pelvic Harness, 2 - A:10.5" x B:9.75" (89108-2) (K0108)	\$290

#### Footplates



Formed Footplates, Small, Size 1 and 2 (Pair) (86147)	\$195
Formed Footplates, Medium, Size 3 and 4 (Pair) (86148)	\$195

## Headrest and Mounting Hardware

### Anatomic Headrest

#### Anatomic Head Support

Anatomic Head Support, Small (5099001) (E0955)	\$302
Anatomic Head Support, Medium (5099002) (E0955)	\$302
Anatomic Head Support, Large (5099003) (E0955)	\$302

#### Head Support Cover

Anatomic Head Support Cover, Small (9923-83)	N/C
Anatomic Head Support Cover, Medium (9945-83)	N/C
Anatomic Head Support Cover, Large (9967-83)	N/C

#### Headrest Mounting System

*Required selection. Required to mount headrest.*

Headrest Mounting Hardware Multi-Adjustable, Short (98619-25) (E1028)	\$245
Headrest Mounting Hardware Multi-Adjustable, Long (98619-26) (E1028)	\$245

### Adjustable Lateral Headrests

#### Adjustable Lateral Headrest

Adjustable Width Headrest, Large Laterals, Black Cover, Short Mounting Hardware (9973792-85)	\$325
Adjustable Width Headrest, Large Laterals, Black Cover, Long Mounting Hardware (9973793-85)	\$245
Adjustable Width Headrest, Large Laterals, with Speakers, Black Cover, Short Mounting Hardware (9971792-85)	\$431
Adjustable Width Headrest, Large Laterals, with Speakers, Black Cover, Long Mounting Hardware (9971793-85)	\$431
Adjustable Width Headrest, Four Laterals, Short Mounting Hardware (9975792-85)	\$403
Adjustable Width Headrest, Four Laterals, Long Mounting Hardware (9975793-85)	\$403
Adjustable Width Headrest, Two Laterals, Short Mounting Hardware (9977792-85)	\$403
Adjustable Width Headrest, Two Laterals, Long Mounting Hardware (9977793-85)	\$403

### Mounting Hardware

Small Mounting Bar	\$245
Long Mounting Bar	\$245

### Headrest Accessories

#### Headband System

Headband System for Anatomic Head Support Only (99500) (E0955)	\$273
--	-------



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



### Headrest Mount

Universal Headrest Mount (87590) (E1028) \_\_\_\_\_ \$101

## Accessories

### Trays

#### Two Rod Tray System

Tray, Two Rod System, Grey, Size 1 and 2, (Each) (8910725-1) (E0950) \$136  
 Tray, Two Rod System, Grey, Size 2 and 3, (Each) (8910725-2) (E0950) \$136  
 Tray, Two Rod System, Grey, Size 4, (Each) (8910725-3) (E0950) \_\_\_\_\_ \$136

#### Two Rod Tray System Hardware

Swing Away Removable Hardware, Size 1 and 2 (Each) (X8910725-12) (E1028) \_\_\_\_\_ \$225  
 Swing Away Removable Hardware, Size 2 and 3 (Each) (X8910725-22) (E1028) \_\_\_\_\_ \$225  
 Swing Away Removable Hardware, Size 4 (Each) (X8910725-32) (E1028) \_\_\_\_\_ \$225

### Tray Accessories

*Not available with Size 4*

#### Toy Bar

Toy Bar for Tray (8910740) \_\_\_\_\_ \$267

#### Arm Pads, Mounts on Top of Tray

Arm Pads, Mounts on Top of Tray, Size 1, 2, 3 (8910730) \_\_\_\_\_ \$343

### Canopy and Hardware

#### Canopy

Canopy, Black, Size 1 and 2 (8740110) \_\_\_\_\_ \$390  
 Canopy, Black, Size 3 and 4 (8740120) \_\_\_\_\_ \$390

#### Canopy Hardware

Holders for Canopy (9589100) \_\_\_\_\_ \$55

## Footrests for Mobility Bases

*Footrest attaches to the mobility bases. If mobility base is desired with a footrest make a selection here, then select desired mobility base below.*

### Footrests and Leg and Foot Supports

*Select footrest type and desired mounting bar*

#### Footrest Type

Angle and Height Adjustable Footplate, Small, 12" x 7.5" (95145-1) (K0040) \_\_\_\_\_ \$281  
 Angle and Height Adjustable Footplate, Medium, 14" x 9" (95145-2)(K0040) \_\_\_\_\_ \$281  
 Angle and Height Adjustable Footplate, Large, 16" x 9.5" ( 95145-3) (K0040) \_\_\_\_\_ \$281

#### Footrest - Multi-Adjustable

*Mounting hardware included. Do not quote mounting bar. Not available with a Foot Box.*

Multi-Adjust Foot Support, Small - 6" - 10.5" (8910852-1) \_\_\_\_\_ \$431  
 Multi-Adjust Foot Support, Medium - 6" - 14.5" (8910852-2) \_\_\_\_\_ \$431  
 Multi-Adjust Foot Support, Large - 6" - 17.5" (8910852-3) \_\_\_\_\_ \$431

### Foot Box Kits

*Select desired size for Footbox and corresponding Foot Plate Padding, Strap for Foot Box is required for all sizes. not available with Multi-Adjustable Footrests.*

#### Foot Box Kit

Foot Box Kit, Size 1 Plate (945145-1 BOX) \_\_\_\_\_ \$240  
 Foot Box Kit, Size 2 Plate (945145-2 BOX) \_\_\_\_\_ \$255  
 Foot Box Kit, Size 3 Plate (945145-3 BOX) \_\_\_\_\_ \$275

#### Strap for foot box

Strap for Foot Box (95145-20) \_\_\_\_\_ \$98

#### Foot Plate Padding

Foot Plate Padding, Size 1 (95145-1PAD) \_\_\_\_\_ \$120  
 Foot Plate Padding, Size 2 (95145-2PAD) \_\_\_\_\_ \$165  
 Foot Plate Padding, Size 3 (95145-3PAD) \_\_\_\_\_ \$205

#### Padding for foot plate rod

Padding for Foot Plate Rod (95145-150) \_\_\_\_\_ \$64

#### Mounting Bar

*Required for Angle and Height Adjustable footplate*

Short Mounting Bar - 6-10.5" (8910821-1) \_\_\_\_\_ N/C  
 Medium Mounting Bar - 6-14.5" (8910821-2) \_\_\_\_\_ N/C  
 Long Mounting Bar - 6-17.5" (8910821-3) \_\_\_\_\_ N/C

#### Width Extensions

*Available only with Multi-Adjustable Foot Supports*

Width Extension - Multi Adjust \_\_\_\_\_ \$82

### Footrest Accessories

Neoprene Sandals - Small (85696-1) \_\_\_\_\_ \$138  
 Neoprene Sandals - Medium (85696-2) \_\_\_\_\_ \$138  
 Formed Shoeholder - Small (86147) \_\_\_\_\_ \$176  
 Formed Shoeholder - Medium (86148) \_\_\_\_\_ \$176  
 Footstraps w/ fix lock (85096) (E0951) \_\_\_\_\_ \$39  
 Ankle Straps, X-Small - C: 6.25"-9" (85745-XS) \_\_\_\_\_ \$197  
 Ankle Straps, Medium - C: 9"-12.25" (85745-M) \_\_\_\_\_ \$197  
 Ankle Straps, X-Large - C: 11.5"-15" (85745-XL) \_\_\_\_\_ \$197



*Attaches to mounting bar*

Calf Support for Centrally Mounted Foot Plate (85131) \_\_\_\_\_ \$245



Phone (800) 981-1540 | Fax (715) 254-0996  
 sales@kimobility.com | www.kimobility.com



## High-low:xo Indoor Base

### High-low:xo Indoor Bases

#### Manual Lifts

High-Low:x Manual Lift, Indoor, Size 1 (89010522) (K0108)	\$2153
High-Low:x, Manual Lift, Indoor, Size 2 (89020522) (K0108)	\$2310
High-Low:x, Hydraulic Lift, Indoor, Size 3 (89030542)(K0108)	\$3468

#### Power Lifts

High-Low:x Indoor, Power Lift, Size 3 (89030532)(K0108)	\$4043
---	--------

### High-low:xo Indoor Bases Accessories

#### Adaptor Mounting Plate

Adaptor Plate for x:panda sz 4 on High-Low:xo 3 (8911610-99)	\$207
--	-------

#### Anti-Tip

Anti-Tip for Indoor High-Low:x Indoor (890007)	\$243
--	-------

#### Push Brace

*Required selection for indoor High:low Base*

Push Brace, Fixed Angle (890009-04)	\$190
Push Brace, Angle Adjustable (8910555)	\$456

## Chilla Mobility Base

### Chilla Frame

#### Chilla Folding Frame

*Only available with size 1 and size 2*

Chilla, Size 2 Frame Solid Adjustable Suspension (881121110) E1234	\$2630
--	--------

#### Shape Adapter

x:panda Shape adapter Plate (881091) K0208	\$875
--	-------

### Chilla Frame Accessories

Chilla Transit Fitting for Shape Adapter (881089)	\$224	Vent Tray (881095)	\$330
IV Pole (881097)	\$475	02 Holder (881096)	\$175

## Multi-Frame Mobility Base

### R82 Multi Frame

#### Multi Frame Configuration

Multi Frame Size 3, Silver, 12.5" Solid Tire, Std. Wheel Locks (973140210)	\$2778
Multi Frame Size 3, Silver, 12.5" Solid Tire, Drum Break (973141210)	\$2778
Drum Break for Multi Frame	\$326

### Multi-Frame Accessories

#### Adapter

Adaptor for x:panda on Multi Frame, Size 1-3 (8911610-98)	\$208
Adaptor for x:panda on Multi Frame, Size 4 (8911610-99)	\$207

#### Accessories

Anti-Tips (970957) (E0971)	\$172
Transit (970905)	\$237
Vent Tray for Multi Frame (970815) (E1029)	\$532



Phone (800) 981-1540 | Fax (715) 254-0996  
sales@kimobility.com | www.kimobility.com

x:panda

## Ki Mobility Bases

*Separate Order form for base must be filled out. Order form is located at Kimobility.com. Footrests are not provided with the base. If footrest are desired, selections must be made above in the "Mobility Base Footrest" section.*

### Ki Mobility Bases

#### **Tyke x:panda**

*Only available with Size 1 and Size 2. Prices will be depicted on the Kimobility.com order form. Must be configured separately.*

Tyke for x:panda \_\_\_\_\_ N/C

#### **Little Wave Flip for x:panda**

*Prices will be depicted on the Kimobility.com order form. Must be configured separately.*

Flip for x:panda \_\_\_\_\_ N/C