



Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

**x:panda shape**  
advanced

Company: \_\_\_\_\_ Ship To: \_\_\_\_\_  
 Account Number: \_\_\_\_\_ Date: \_\_\_\_\_ Address: \_\_\_\_\_  
 Quote: PO Number: \_\_\_\_\_ City: \_\_\_\_\_ State: \_\_\_\_\_  
 ATP: \_\_\_\_\_ Zip: \_\_\_\_\_ Phone: \_\_\_\_\_  
 Marked For: \_\_\_\_\_ Email/Fax: \_\_\_\_\_  
 Associated Quote/Order#: \_\_\_\_\_ Effective: 2/3/2025

**x:panda Shape - Advanced**

	Infant	1	2
Effective seat width	8.5" 220 mm	8.5"/10.25"/11.75" 220/260/300 mm	11"/12.5"/14.25"/15.75" 280/320/360/400 mm
Effective seat depth	4.75-6.25" 120-160 mm	6.25-12.5" 160-320 mm	9.5-16.5" 240-420 mm
Backrest height	10.75-14.25" 270-360 mm	11-14.5"/12.5-16.5"/14.25-18.5" 280-370/320-420/360-470mm	13.5-17.25"/15-19.25"/17-22" 340-440/380-490/430-560 mm
Max user mass	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
Max user mass in transportation	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
Backrest angle	-10° - 30°	-10° - 30°	-10° - 30°
Back width, top of backrest	7.5" 190 mm	7"/8.5"/10.25"/11.75"/13.5" 180/220/260/300/340 mm	7"/8.5"/10.25"/11.75"/13.5"/15" 180/220/260/300/340/380 mm
Overall width	9.5" 240 mm	11.75"/15.75"/19.75" 300/400/500 mm	13.75"/18"/22"/26" 350/460/560/660 mm
Mass device	13.2 lb 6 kg	14.3 lb 6.5 kg	24.3 lb 11 kg

Measurements taken with both legs in max abduction. Size 1-3 measurements are based on Advanced Seat.

**Shape Base Options**

Please see measurement guide for growth information. A selection must be made in each section below.

**Seat Base**

Includes base and plate

**Seat Base Size**

Size 1 Advanced Seat (E2292) \_\_\_\_\_ \$1122  
 Size 2 Advanced Seat (E2292) \_\_\_\_\_ \$1122

**Recline Options**

Selection required

Reclining Dynamic Back (K0108) \_\_\_\_\_ \$709  
 Reclining Adjustable Dynamic Back (K0108) \_\_\_\_\_ \$913

HCPCS Codes provided are based on code verification by the PDAC or our interpretation of the code definitions.  
 Proper use of HCPCS Codes does not ensure coverage or payment. For coverage information, verify the appropriate payer's coverage policy.  
 DCN0626.4 Pricing and specifications subject to change without notice.



Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

x:panda shape  
advanced

## Back and Cushions

### Upper Back Size

*"Select Upper Back Base then select Upper and Lower Back Cushion and Covers. Includes quick release detachment (#9050374)*



#### Size 1 Base

Size 1 Base Upper - H: 7" (905073) (E2291)	\$255
<i>Total back height range 12.25" - 16.25"</i>	
Size 1 Base Upper - H: 8.75" (905074) (E2291)	\$255
<i>Total back height range 13.75" - 18.50"</i>	

#### Size 2 Base

Size 2 Base Upper - H: 7" (905073) (E2291)	\$255
<i>Total back height range 13.5" - 17.25"</i>	
Size 2 Base Upper - H: 8.75" (905074) (E2291)	\$255
<i>Total back height range 15" - 19.75"</i>	
Size 2 Base Upper - H: 10.75" (905075) (E2291)	\$255
<i>Total back height range 17" - 21.75"</i>	

### Upper Back Cushion and Cover

*Upper back cushion height must match Upper Back Base height*

#### Upper Back Cushion

Upper Back Cushion - W: 7" x H: 7" (902103)	N/C
Upper Back Cushion - W: 8.75" x H: 7" (902104)	N/C
Upper Back Cushion - W: 10.25" x H: 7" (902105)	N/C
Upper Back Cushion - W: 8.75" x H: 8.75" (902106)	N/C
Upper Back Cushion - W: 10.25" x H: 8.75" (902107)	N/C
Upper Back Cushion - W: 11.75" x H: 8.75" (902108)	N/C
Upper Back Cushion - W: 13.5" x H: 8.75" (902109)	N/C
Upper Back Cushion - W: 10.25" x H: 10.75" (902110)	N/C
Upper Back Cushion - W: 11.75" x H: 10.75" (902111)	N/C
Upper Back Cushion - W: 13.5" x H: 10.75" (902112)	N/C
Upper Back Cushion - W: 15" x H: 10.75" (902113)	N/C

#### Upper Back Cover

*Selection required*

Red (-50)	N/C
Grey (-60)	N/C
Black Wipeable (-70)	\$21

### Lower Back Cushion and Cover

*The width measurements of the cushion are not actual measurements of the cushion. The number matches the lower back of the seating system. If the seating system is configured without hip and pelvis support, consider going up a size in width when choosing the lower back cushion.*

#### Lower Back Cushion

Size 1 Lower Back Cushion - W: 8.75" x H: 4.75" (902120)	N/C
Size 1 Lower Back Cushion - W: 10.25" x H: 4.75" (902121)	N/C
Size 1 Lower Back Cushion - W: 11.75" x H: 4.75" (902122)	N/C
Size 2 Lower Back Cushion - W: 11" x H: 5.5" (902123)	N/C
Size 2 Lower Back Cushion - W: 12.5" x H: 5.5" (902124)	N/C
Size 2 Lower Back Cushion - W: 14.25" x H: 5.5" (902125)	N/C
Size 2 Lower Back Cushion - W: 15.75" x H: 5.5" (902126)	N/C

#### Lower Back Cover

*Selection Required*

Red (-50)	N/C
Grey (-60)	N/C
Black Wipeable (-70)	\$21

### Daily Adjustment Kit for Upper Back

*Allows for tool-less adjustment of back height. Ideal for frequent adjustment.*



Daily Adjustment Kit - Upper Back (905081)	\$21
--	------

## Spacer Blocks

### Spacer Block Holder without Cushion

*Cushion width must match spacer block holder width*

Spacer Block Holder, Size 1 - Back W: 7", 8.75", 10.25"	\$179
Spacer Block Holder, Size 2 - Back W: 8.75", 10.25", 11.75", 13.5"	\$179

### Spacer Block Cushion Cover

Black (-80)	N/C
Black Wipeable (-70)	N/C

### Spacer Block Cushion

Spacer Block Cushion Size 1 - W: 7" x H: 2.25" (903011)	N/C
Spacer Block Cushion Size 2 - W: 8.75" x H: 2.25" (903012)	N/C
Spacer Block Cushion Size 3 - W: 10.25" x H: 2.25" (903013)	N/C
Spacer Block Cushion Size 4 - W: 11.75" x H: 2.25" (903014)	N/C
Spacer Block Cushion Size 5 - W: 13.5" x H: 2.25" (903015)	N/C





Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

x:panda shape  
advanced

## Cushions and Covers

### Seat Cushion

#### Seat Cushion - Split

*Required when using Quick Release Pommel*

Split Cushion, Size 1- W: 8.75" x L: 19.75" (901113)	N/C
Split Cushion, Size 1- W: 10.25" x L: 19.75" (901114)	N/C
Split Cushion, Size 1- W: 11.75" x L: 19.75" (901115)	N/C
Split Cushion, Size 2- W: 12.5" x L: 25.5" (901117)	N/C
Split Cushion, Size 2- W: 14.25" x L: 25.5" (901118)	N/C
Split Cushion, Size 2- W: 15.75" x L: 25.5" (901119)	N/C

#### Color

*Required selection for cover color*

Red (-50)	N/C
Grey (-60)	N/C
Black Wipeable (-70)	\$21

## Headrest and Hardware Options

### Headrest Style

*Includes Cushion*

#### Type of Head Support

Curved, Small - Black Wipeable - W: 7", H: 5" (E0955)	\$220
Curved, Small - Black - W: 7", H: 5" (E0955)	\$220
Curved, Medium - Black Wipeable (E0955)	\$240
Curved, Medium - Black (E0955)	\$240
Anatomic, Small - Black Wipeable (E0955)	\$271
Anatomic, Small - Black (E0955)	\$271
Anatomic, Medium - Black Wipeable (E0955)	\$281
Anatomic, Medium - Black (E0955)	\$281
Anatomic, Large - Black Wipeable (E0955)	\$291
Anatomic, Large - Black (E0955)	\$291
Adjustable with Large Laterals - Black (E0955)	\$325
Adjustable with Large Laterals and Speakers - Black (E0955)	\$431
Adjustable with Four Laterals - Black (E0955)	\$403

### Mounting Hardware for Curved and Anatomic Headrest

*Required with Curved and Anatomic Head support*

#### Headrest Mounting Hardware

Mounting Brackets, Standard (E1028)	\$245
5.5" of forward/backward adjustment	
Mounting Brackets, Extended (E1028)	\$255
8.75" of forward/backward adjustment	
Mounting Brackets, 3D Adjustment (E1028)	\$266
8.75" of forward/backward adjustment	

### Mounting Hardware for Adjustable Headrest Style

*Required with Adjustable Head Support*

Short Mounting Hardware (E1028)	\$245
Suitable for upper back heights of 7", 8.75", 10.75" and 13.5"	
Long Mounting Hardware (E1028)	\$245
Suitable for upper back heights of 10.75" and 13.5"	

### Headrest Options



#### Daily Adjustment Kits

*Creates toolless adjustability, For angle adjustment of head support. For use with Anatomic and Curved Headrests only.*

Daily Adjustment Kit Standard - Set of Two	\$31
Daily Adjustment Kit Extended - Set of Three	\$46

### Additional Headrest Options

Headband System	\$245
For use with Anatomic Headrest	
Universal Headrest Mount	\$66
Required if using a third party headrest	

## Arm Supports

### Arm Supports

#### Arm Supports Size

*Cover not included, select separately*

Size 2 - W: 3" L: 10.25" H: 4" - 6.5" (8910616-1) (E0973)	\$332
Size 3 - W: 3" L: 12.25" H: 4" - 8.75" (8910617-1) (E0973)	\$332

#### Arm Support Cover

Black	\$77
Black Wipeable	\$77

#### Daily Adjustment Kit

*For toolless height adjustment of arm supports*

Height Adjustment Kit of Arm Supports (8910615-11)	\$21
--	------



#### Upper Arm Support - Fixed

*Includes Black Wipeable cover. Not compatible with Swing-Away Lateral Supports.*

Size 1 - L: 5.5" H: 7" (905017-70) (K0108)	\$255
For size 1 Base	
Size 2 - L: 7.5" H: 10.75" (905018-70) (K0108)	\$255
For size 2 Base	

#### Upper Arm Supports - Adjustable

*Not compatible with Swing-Away Lateral Supports.*

Size 1 - L: 5.5" H: 7" (905019-70) (K0108)	\$342
For Size 1 Base	
Size 2 - L: 7.5" H: 10.75" (905020-70) (K0108)	\$342
For size 2 Base	



Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

x:panda shape  
advanced

## Knee and Postural Supports

### Abduction and Adduction

#### Inner Knee Support



For Size 1 and 2 Black - L: 7.5", H: 7.5" (905001-80) (E0957)	\$296
For Size 1 and 2 Black Wipeable - L: 7.5", H: 7.5" (905001-70) (E0957)	\$296
For Size 2 Black - L: 9.5", H: 9.5" (905002-80) (E0957)	\$296
For Size 2 Black Wipeable - L: 9.5", H: 9.5" (905002-70) (E0957)	\$296

#### Outer Knee Support



Standard Set for Size 1, Black - L: 6", H: 3.25" (905008-80) (E0956)	\$215
Standard Set for Size 1, Black Wipeable - L: 6", H: 3.25" (905008-70) (E0956)	\$215
Standard Set for Size 2, Black - L: 7.75", H: 3.5" (905009-80) (E0956)	\$215
Standard Set for Size 2, Black Wipeable - L: 7.75", H: 3.5" (905009-70) (E0956)	\$215
Flip Down Set for Size 2, Black - L: 7.75", H: 7.75" (905010-80) (E0956)	\$255
Flip Down Set for Size 2, Black Wipeable - L: 7.75", H: 7.75" (905010-70) (E0956)	\$255

#### Quick Release Pommel

For Size 1 Small, Black - W: 2.25", L: 6.25", H: 7.5" (905077-80) (E0957)	\$357
For Size 1 Small, Black Wipeable - W: 2.25", L: 6.25", H: 7.5" (905077-70) (E0957)	\$357
For Size 2 Small, Black - W: 2.25", L: 6.25", H: 7.5" (905078-80) (E0957)	\$357
For Size 2 Small, Black Wipeable - W: 2.25", L: 6.25", H: 7.5" (905078-70) (E0957)	\$357
For Size 2 Large, Black - W: 3.25", L: 7.5", H: 9" (905082-80) (E0957)	\$357
For Size 2 Large, Black Wipeable - W: 3.25", L: 7.5", H: 9" (905082-70) (E0957)	\$357

### Postural Supports

#### Width Adjustable Pelvis Supports

*Mounted from the back*



Width Adjustable Pelvis Supports, Small Black - L: 4", H: 4.25" (905041-80) (E0956)	\$240
Width Adjustable Pelvis Supports, Small Black Wipeable - L: 4", H: 4.25" (905041-70) (E0956)	\$240
Width Adjustable Pelvis Supports, Large Black - L: 5.5", H: 5" (905042-80) (E0956)	\$240
Width Adjustable Pelvis Supports, Large Black Wipeable - L: 5.5", H: 5" (905042-70) (E0956)	\$240

#### Width Adjustable Hip and Pelvis Supports

*Mounted from the sides*



Width Adjustable Hip Supports, Small Black - L: 4" H: 4.25" (905044-80) (E0956)	\$342
Width Adjustable Hip Supports, Small Black Wipeable - L: 4" H: 4.25" (905044-70) (E0956)	\$342
Width Adjustable Hip Supports, Large Black - L: 5.5" H: 5" (905045-80) (E0956)	\$342
Width Adjustable Hip Supports, Large Black Wipeable - L: 5.5", H: 5" (905045-70) (E0956)	\$342

#### Adjustable Anti-Thrust Wedges

Short - W: 4.5", L: 2.75-4.25", H: .75" (905052) (K0108)	\$113
<i>For Size 1</i>	
Medium - W: 5.25", L: 4.25-6.75", H: 1.25" (905053) (K0108)	\$113
<i>For Size 1</i>	
Long - W: 5.25", L: 4.75-7.75", H: 1.75" (905054) (K0108)	\$113
<i>For Size 1</i>	
Medium - W: 6", L: 5-8.75", H: 1.25" (905055) (K0108)	\$113
<i>For Size 2</i>	
Long - W: 6", L: 5.50-9.5", H: 2" (905056) (K0108)	\$113
<i>For Size 2</i>	

## Trays and Accessories



Shape 5 Tray with Edge Size 1 (8910725-1) (E0950)	\$327
<i>For Size 1. Fits Size 2 widths 11", 12.5", 14.25". Double Rod, Grey ABS.</i>	



Shape 6 Tray with Edge Size 2 (8910725-2) (E0950)	\$327
<i>For Size 2. Double rod, Grey ABS.</i>	



Shape 8 Tray with Edge Size 1 (8976350) (E0950)	\$327
---	-------

Padding for Size 1 Tray, Shape 5 Plate, Over and Under (8910725-10) (E1028)	\$434
Padding for Size 1 Tray, Top, Shape 5 Plate, Over (8910725-10T) (E1028)	\$306
Padding for Size 2 Tray, Shape 6 Plate, Over and Under (8910725-20) (E1028)	\$434
Padding for Size 2 Tray, Top, Shape 6 Plate, Over (8910725-20T) (E1028)	\$306
Arm Pads for Tray (8910730)	\$309
<i>For shape 5 and 6</i>	
Toy Bar (8910740)	\$245



Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

# x:panda shape advanced

## Accessory Holders and Brackets

Remote Control Holder for Size 1 and 2 - Arm Support Mounted (905030-1)	\$67
Remote Control Holder for Size 3 - Arm Support Mounted (905030-2)	\$67
Remote Control Holder - Back Mounted (905030)	\$128

## Mounting Hardware

Bracket for Communication Aid, Size 1 and 2 - Arm Support Mounted (905033-1)	\$80
Bracket for Communication Aid, Size 3 - Arm Support Mounted (905033-2)	\$80

## Canopy

*Canopy does not come with the brackets. Both the Canopy and brackets must be quoted when mounting to shape.*



### Canopy Without Brackets

Black Canopy - Size 1 (8740110)	\$332
Black Canopy - Size 2 (8740120)	\$332

## Canopy Brackets

Canopy Brackets with Push Handle for Size 1 (905039-1)	\$199
Canopy Brackets with Push Handle for Size 2 (905039-2)	\$199

## Vest and Belts

*More options can be found at etac.com*

### Vests

#### Cross Vests



Cross Vest X-Small - S-P: 7.75" - 12.5" (85720-XS) (E0960)	\$202
Cross Vest Small - S-P: 10.25" - 13.5" (85720-S) (E0960)	\$202
Cross Vest Medium - S-P: 13.5" - 17.75" (85720-M) (E0960)	\$202
Cross Vest Large - S-P: 17.75" - 20.75" (85720-L) (E0960)	\$202
Cross Vest X-Large - S-P: 20.75" - 24.5" (85720-XL) (E0960)	\$202



#### Torso Vests



Torso Vest Small - S-P: 10.75" - 12.5" (85705-S) (E0960)	\$214
Torso Vest Medium - S-P: 12.5" - 16" (85705-M) (E0960)	\$214
Torso Vest Large - S-P: 16" - 18.5" (85705-L) (E0960)	\$214
Torso Vest X-Large - S-P: 18.5" - 20.5" (85705-XL) (E0960)	\$214

#### Fixlocks Bracket



Higher Brackets (8910260)	\$77
<i>Bracket for higher strap guidance. Includes set of fixlocks. Requires the wider brackets as well.</i>	
Wider Brackets (905068)	\$57

### Pelvic Belts and Harnesses

#### 4-Point Belts



4-Point Padded Hip Belts 1" X-Small - H: 6.75" - 12.5" (85725-XS) (E0960)	\$143
4-Point Padded Hip Belts 1" Small - H-H: 8.25" - 13.75" (85725-S) (E0960)	\$143
4-Point Padded Hip Belts 1" Medium - H: 9.75" - 16.25" (85725-M) (E0960)	\$143
4-Point Padded Hip Belts 1" Large - H-H: 12.5" - 17.25" (85725-L) (E0960)	\$143



#### Pelvic Harness

Pelvic Harness, 0 - A:6.5" x B:7.5" (89108-0)	\$261
Pelvic Harness, 1 - A:8.5" x B:8.75" (89108-1)	\$261
Pelvic Harness, 2 - A:10.5" x B:9.75" (89108-2)	\$261



#### 2-Point Pelvic Belt

2-Point Pelvic Belt Black - L: 13-14.5" (85424-1) (E0978)	\$104
---	-------



#### Pelvis Belt

Pelvis Belt, X-Small - H-H: 6.25" - 7.75" (85735-XS) (E0978)	\$261
Pelvis Belt, Small - H-H: 7.75" - 9.5" (85735-S) (E0978)	\$261
Pelvis Belt, Medium - H-H: 9.5" - 11" (85735-M) (E0978)	\$261



Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

x:panda shape  
advanced

## Lateral Support Cushion and Fittings

### Lateral Supports

If choosing Laterals, both the Left and Right cushions must be selected.

#### Left Cushion

Left Rectangular Small - L: 4.25" x H: 3.25" (E0956)	\$113
Left Rectangular Medium - L: 6" x H: 3.5" (E0956)	\$113
Left Rectangular Large - L: 7.5" x H: 3.5" (E0956)	\$113

#### Right Cushion

Right Rectangular Small - L: 4.25" x H: 3.25" (E0956)	\$113
Right Rectangular Medium - L: 6" x H: 3.5" (E0956)	\$113
Right Rectangular Large - L: 7.5" x H: 3.5" (E0956)	\$113

#### Type of Fitting

Selection required, fittings are sold as a pair.

Fixed Fitting (E1028)	\$169
Swing-Away and Angle Adjustable Fitting (E1028)	\$255

## Footrests for Mobility Bases

Footrest attached the mobility bases. If mobility base are desired with a footrest make a selection here and select desired mobility base below.

### Footrests

Select footrest type and desired mounting bar

#### Footrest Type

Angle and Height Adjustable Footplate, Small, 12" x 7.5" (95145-1) (K0040)	\$281
Angle and Height Adjustable Footplate, Medium, 14" x 9" (95145-2)(K0040)	\$281
Angle and Height Adjustable Footplate, Large, 16" x 9.5" (95145-3) (K0040)	\$281

#### Footrest - Multi-Adjustable

Mounting hardware included. Do not quote mounting bar. Not available with Foot Box,

Multi-Adjust Foot Support, Small - 6" - 10.5" (8910852-1)	\$431
Multi-Adjust Foot Support, Medium - 6" - 14.5" (8910852-2)	\$431
Multi-Adjust Foot Support, Large - 6" - 17.5" (8910852-3)	\$431

### Foot Box Kits

Select desired size for Foot Box and corresponding Foot Plate Padding. Strap for Foot Box is required with all sizes. Not available with Multi Adjustable Footrests.

#### Foot Box Kit

Foot Box Kit, Size 1 Plate (945145-1 BOX)	\$240
Foot Box Kit, Size 2 Plate (945145-2 BOX)	\$255
Foot Box Kit, Size 3 Plate (945145-3 BOX)	\$275

#### Strap for foot box

Strap for Foot Box (95145-20)	\$98
-------------------------------	------

#### Foot Plate Padding

Foot Plate Padding, Size 1 (95145-1PAD)	\$120
Foot Plate Padding, Size 2 (95145-2PAD)	\$165
Foot Plate Padding, Size 3 (95145-3PAD)	\$205

#### Padding for foot plate rod

Padding for Foot Plate Rod (95145-150)	\$64
--	------

#### Mounting Bar

Required for Angle and Height Adjustable Footplates

Short Mounting Bar - 6-10.5" (8910821-1)	N/C
Medium Mounting Bar - 6-14.5" (8910821-2)	N/C
Long Mounting Bar - 6-17.5" (8910821-3)	N/C

#### Width Extensions

Available only with Multi-Adjustable Foot Supports

Width Extension - Multi Adjust	\$82
--------------------------------	------

### Footrest Accessories

Neoprene Sandals - Small (85696-1)	\$138
Neoprene Sandals - Medium (85696-2)	\$138
Formed Shoeholder - Small (86147)	\$176
Formed Shoeholder - Medium (86148)	\$176
Footstraps w/ fix lock (85096) (E0951)	\$39
Ankle Straps, X-Small - C: 6.25"-9" (85745-XS)	\$197
Ankle Straps, Medium - C: 9"-12.25" (85745-M)	\$197
Ankle Straps, X-Large - C: 11.5"-15" (85745-XL)	\$197



Attaches to mounting bar

Calf Support for Centrally Mounted Foot Plate (85131)	\$245
---	-------





Phone (877) 227-0097 | Fax (704) 882-0751  
sales.us@etac.com

x:panda shape  
advanced

## Chilla Mobility Base

### Chilla Frame

#### Chilla Folding Frame

Chilla, Size 2 Frame Solid Adjustable Suspension (881121110) E1234 \_\_\_\_\_ \$2630

#### Shape Adapter

x:panda Shape adapter Plate (881091) K0208 \_\_\_\_\_ \$875

### Chilla Frame Accessories

Chilla Transit Fitting for Shape Adapter (881089) \_\_\_\_\_ \$224      O2 Holder (881096) \_\_\_\_\_ \$175

Vent Tray (881095) \_\_\_\_\_ \$330      IV Pole (881097) \_\_\_\_\_ \$475

## High-low:xo Mobility Base and Options

### High-low:xo

#### Manual Lifts

Manual Lift - indoor Size 1 (K0108) \_\_\_\_\_ \$1937

Manual Lift - indoor, Size 2 (K0108) \_\_\_\_\_ \$2079

#### Outdoor High-low:xo

Outdoor- manual lift, Size 2 (89120152) \_\_\_\_\_ \$3264

Outdoor- power lift, Size 2 \_\_\_\_\_ \$4182

#### Power Lifts

Indoor- power lift, Size 1 (K0108) \_\_\_\_\_ \$3113

Indoor- power lift, Indoor, Size 2 (K0108) \_\_\_\_\_ \$3342

#### Anti-Tips on the High-low:xo

Anti-tips - Indoor (890007) \_\_\_\_\_ \$219

Anti-tips - Outdoor (890008) \_\_\_\_\_ \$235

#### Transit - Outdoor High-low:xo

Transit (890005) \_\_\_\_\_ \$255

### Push Brace for High-low:xo - Indoor

#### Required selection

Fixed Angle Push Brace (890009-04) \_\_\_\_\_ \$194

Angle Adjustable Push Brace (8910555) \_\_\_\_\_ \$411

## Multi Frame Mobility Base and Options

### Multi Frame Base Options

Multi Frame Size 3, Silver, 7" x 2", 12.5" Solid Tire, Standard Wheel \_\_\_\_\_ \$2501

Locks \_\_\_\_\_ \$2796

Multi Frame Size 3, Silver, 7" x 2", 12.5" Solid Tire, Drum Brake \_\_\_\_\_ \$2796

#### Anti- Tip Option

Anti-tips \_\_\_\_\_ \$155

#### Transit Option

Transit for Multi Frame \_\_\_\_\_ \$214

### Multi-Frame Accessories

Adaptor for x:panda on Multi Frame - for x:panda sz 1 & 2 \_\_\_\_\_ \$187      Vent Tray \_\_\_\_\_ \$490

O2 Holder \_\_\_\_\_ \$373

## Ki Mobility Bases and Options

### Ki Mobility Bases

*Separate Order form for base must be filled out. Order form is located at Kimobility.com. Footrests are not provided with the base. If footrest are desired, selections must be made above in the "Mobility Base Footrest" section.*

#### Tyke x:panda

*Prices will be depicted on the Kimobility.com order form. Must be configured separately.*

Tyke for x:panda \_\_\_\_\_ N/C

#### Little Wave Flip for x:panda

*Prices will be depicted on the Kimobility.com order form. Must be configured separately.*

Flip for x:panda \_\_\_\_\_ N/C