



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

Company: _____ Ship To: _____
 Account Number: _____ Date: _____ Address: _____
 Quote: PO Number: _____ City: _____ State: _____
 ATP: _____ Zip: _____ Phone: _____
 Marked For: _____ Email/Fax: _____
 Associated Quote/Order#: _____ Effective: 2/3/2025

	Infant	1	2
Effective seat width	8.5" 220 mm	8.5"/10.25"/11.75" 220/260/300 mm	11"/12.5"/14.25"/15.75" 280/320/360/400 mm
Effective seat depth	4.75-6.25" 120-160 mm	6.25-12.5" 160-320 mm	9.5-16.5" 240-420 mm
Backrest height	10.75-14.25" 270-360 mm	11-14.5"/12.5-16.5"/14.25-18.5" 280-370/320-420/360-470mm	13.5-17.25"/15-19.25"/17-22" 340-440/380-490/430-560 mm
Max user mass	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
Max user mass in transportation	77.2 lb 35 kg	77.2 lb 35 kg	143.3 lb 65 kg
Backrest angle	-10° - 30°	-10° - 30°	-10° - 30°
Back width, top of backrest	7.5" 190 mm	7"/8.5"/10.25"/11.75"/13.5" 180/220/260/300/340 mm	7"/8.5"/10.25"/11.75"/13.5"/15" 180/220/260/300/340/380 mm
Overall width	9.5" 240 mm	11.75"/15.75"/19.75" 300/400/500 mm	13.75"/18"/22"/26" 350/460/560/660 mm
Mass device	13.2 lb 6 kg	14.3 lb 6.5 kg	24.3 lb 11 kg

Measurements taken with both legs in max abduction. Size 1-3 measurements are based on Advanced Seat.

Seat Bases

Please see measurement guide for growth information. A selection must be made in each section below.

Set Base

Includes base and plate

Size 1 Infant Base _____ \$811
See chart above for size details

Recline Options

Selection Required

Adjustable Stepless Recline, Toolless (K0108) _____ \$434
Adjustable Stepless Recline with Dynamic, Lock of Dynamic, Toolless (K0108) _____ \$709

Back

Upper Back with Quick Release Mount _____ \$255

Cushions and Covers

Seat Cushion

Cushion Color

Cushion including seat, back, seat pit, upper arm cushions. Order covers for side support belt separately. Required Selection.

Pink	N/C
Light Blue	N/C
Black Wipeable	\$21

Extra Cushion and Cover

An extra cushion and cover set, in addition to the installed cushion and cover selected above. Comes uninstalled.

Pink	\$332
Light Blue	\$332
Black Wipeable	\$357

Extra Cover only

Select for an extra cover (no cushion). Comes uninstalled. Optional selection.

Pink	\$179
Light Blue	\$179

Headrests and Hardware Options

Headrest Style

Type of Head Support

Includes Short Rod length. Includes Cushion.

Curved, Small - Black Wipeable - W: 7", H: 5" (E0955)	\$220
Curved, Small - Black - W: 7", H: 5" (E0955)	\$220
Anatomic, Small - Black, Wipeable (E0955)	\$271
Anatomic, Small - Black (E0955)	\$271

Headrest Mounting Hardware

2 pivot points, standard 80mm rod, height adjustment rod 905037

Standard Mounting Hardware (E1028)	\$245
------------------------------------	-------



Headrest Options

Knob replaces bolt for tool-less adjustment of headrest angle

Daily Adjustment Kit Standard - Set of Two	\$31
Daily Adjustment Kit Standard - Set of Three	\$46

Vest, Straps and Belts



Side Supports Belts

Infant Side Support Belt, Black Wipeable (E0978)	\$215
Infant Side Support Belt, Black (E0978)	\$215

Torso Vest

Torso Vest X-Small - S-P: 7.75" - 10.75" (E0960)	\$214
Torso Vest Small - S-P: 10.75" - 12.5" (E0960)	\$214

Side Support Extra Cover

Extra Side Support Cover, Black, Wipeable	\$97
Extra Side Support Cover, Black	\$97

Cross Vest

Cross Vest X-Small - S-P: 7.75" - 12.5" (E0960)	\$202
Cross Vest Small - S-P: 10.75" - 13.5" (E0960)	\$202

Padded Hip Belt

4-Point Padded Hip Belt, 1" X-Small - H-H: 6.75" - 12.5" (E0978)	\$143
4-Point Padded Hip Belt, 1" Small - H-H: 8.25" - 13.75" (E0978)	\$143
2-Point Pelvic Belt Black - L: 13-14.5" (E0978)	\$104

Trays and Accessories

Tray and Activities Centres



Tray

Soft tray for Infant Seating, Black - W: 12.5" x L: 7.5"	\$153
--	-------

Canopy

Canopy

Black Canopy - Size 1	\$332
-----------------------	-------

Canopy Brackets

Canopy Brackets with Push Handle for Size 1	\$199
---	-------



Phone (877) 227-0097 | Fax (704) 882-0751
sales.us@etac.com

Footrests for Mobility Bases

Footrest attached the mobility bases. If mobility base are desired with a footrest make a selection here and select desired mobility base below.

Footrests

Footrest Type

Angle and Height Adjustable Footplate, Small, 12" x 7.5"
(95145-1)(K0040) _____ \$281

Mounting Bar

Short Mounting Bar - 6-10.5" (8910821-1) _____ N/C

Footrest Accessories

Neoprene Sandals - Small _____ \$138
Neoprene Sandals - Medium _____ \$138
Footstraps w/ fix lock (E0951) _____ \$39
Ankle Straps, X-Small - C: 6.25" - 9" _____ \$197

Chilla Mobility Base

Chilla Frame

Chilla Folding Frame

Chilla, Size 2 Frame Solid Adjustable Suspension (881121110) E1234 _____ \$2630

Shape Adapter

x:panda Shape adapter Plate (881091) K0208 _____ \$875

Chilla Frame Accessories

Chilla Transit Fitting for Shape Adapter (881089) _____ \$224
Vent Tray (881095) _____ \$330
O2 Holder (881096) _____ \$175
IV Pole (881097) _____ \$475

Mobility Bases

High-low:xo

Manual Lifts

Manual Lift - indoor Size 1 (K0108) _____ \$1937

Anti-Tips on the High-low:xo

Anti-tips - Indoor _____ \$219
Anti-tips - Outdoor _____ \$235

Power Lifts

Indoor- power lift, Size 1 (K0108) _____ \$3113

Push Brace for High-low:xo - Indoor

Required for High-low:xo

Fixed Angle Push Brace _____ \$194
Angle Adjustable Push Brace _____ \$411

Ki Mobility Bases

Separate order form for base must be filled out. Order form is located at www.kimobility.com. Footrests are not provided with base. If footrests are desired selection must be made below.

Tyke for x:panda _____ N/C
Flip for x:panda _____ N/C