



Little Wave Klik

A paediatric wheelchair for our next generation


Little Wave Clik

You are only a “Clik” away from just perfect

Little Wave Clik, the most intuitive paediatric wheelchair, evolves and grows with the child.

Designed for our next generation

An adjustment system designed to provide accurate, easy growth options. A system so easy to adjust with just a “klik” that promotes proper positioning and wheel access so your child gains the active independence and mobility they want. Little Wave Clik has all possibilities to meet the child’s needs.



Evolving and growing with your child.

Seat height, width, depth, Center of Gravity, caster wing position, seat upholstery, seat pans... the most built in adjustment available. With all the possibilities, meeting kids needs has never been so easy.



Expanding your child's world with confidence.

Being able to move in space is an essential phase of neurodevelopment. The challenge is maintaining security and stability. The adjustable dynamic 5th wheel makes this possible, allowing you to adjust the range and spring rate. Whether set up in mid-wheel for the smallest riders or to help promote and teach the balance and wheeling skills that encourage confidence to overcome everyday obstacles and terrain. Be Dynamic.

Features



The adjustment system is designed to provide accurate, easy growth options. Meet kids' needs for independence and mobility with just a "klik".



Whether it's for mom, dad, sister, brother, friend or caregiver the adjustable height Ergo stroller handle is the ergonomic choice for all.




With the accessible space on the frame, you avoid conflict and get access to clamp accessories where you need them.



Compound offset frames ensure your kid has the fit they need. More room for abducted positioning, Ankle Foot Orthoses and braces.

Options



Push Handle
Straight.



Anti Tippers



Footplate
High Mount Angle
Adjustable Flip Under.



Footplate
Angle Adjustable.



Footplate
Angle Adjustable w/
Riser.



Fender
Adjustable aluminium
fender with Foam Pad.



Fender
Adjustable Carbon
fiber fender.



**Personalise with
colour**

Visit etac.com for more options and accessories

Specifications	Imperial	Metric
Transport Weight:	12.5 lbs.*	5.67 kg*
Seat Width:	8-18 in.	20-46 cm
Seat Depth:	8-18 in.	20-46 cm
Front Seat Height:	11.5-20 in.	29-51 cm
Rear Seat Height:	11.5-20.5 in.	29-52 cm
Tubing Size:	1.125 in.	2.85 cm
Weight Capacity:	165 lbs.	75 kg
WC-19 Transit:	Yes	Yes

*12"x12" Frame w/ Adjustable Height Back-12", 4" Poly Casters, Angle Adj. Footplate and Standard Upholstery. Armrest and Rear Wheels Removed.

Founded in 1973, Etac is a world-leading developer and provider of ergonomic assistive devices and patient handling equipment. Our ambition is to provide solutions that optimise quality of life for the individual, family members and caregivers.

For the latest news and continuously updated product information – visit: www.etac.com

Etac

+46 371 587 330
info@etac.se
www.etac.com

etac[®]
Creating Possibilities