

Scandinavian
design with flexible
functionality



Etac Relax shower seat

etac[®]
Creating Possibilities



HEYLAND WHITTLE
Santal de Indes
Santal de Indes
Santal de Indes

etac

etac

Enjoy a relaxing shower

It's lovely having a shower – especially when you can sit down

For most people it's obvious that you should enjoy the time you spend in the shower. For many people, though, standing up in the shower can feel unsafe. However, there may not be space for an ordinary shower stool, which is why we've developed Etac Relax: a neat but secure wall-mounted shower seat. Always close at hand – but because it easily folds up against the wall it doesn't take up unnecessary space. We know how important it is to use contrasting colours to make everything clear in a steam-filled shower. That's why we are now launching Etac Relax in Volcano grey. Perfect for white shower environments.

Attention to detail

Etac Relax shower seat has been developed with great care and attention to detail, in order to make your shower experience more enjoyable.

Scandinavian minimalist design gives you maximum comfort and a more stylish bathroom. Relax has been rigorously tested and approved to meet the stringent requirements of European directives. Design, function and safety are essential for Etac.





Scandinavian design



Ergonomically shaped arm supports with a soft non-slip grip

Etac Relax folds up neatly against the wall



Functional seat shape for stability and easy access

Choose your model

Etac Relax is available in 2 colours and 6 different configurations to meet all needs with regard to making life easier, making a total of 12 versions. Important factors could be how you move in the shower, how much support you need, or the nature of the wall and reducing the load put on it.

12 different versions



Sturdy seat surface

Large seat surface even though the external dimensions are slim.



Non-slip arm supports

The folding arm supports have a non-slip surface, offering a soft and secure grip.




Arm support and back support

A slimline warming back pad with good friction to increase comfort and relaxation.



Versions with supporting legs

The supporting legs are height adjustable to ensure a comfortable seat height and a perfect fit to the floor.



Easy to fit

Quick and easy installation. You only need three screws to fit Etac Relax.

Etac Relax shower seat

Width	40 cm (15¾")
Width with arm supports	53 cm (21")
Depth in up/down position	45 cm (17¾")/8 cm (3")
Maximum user weight	
~without supporting legs	125 kg (276 lbs)
~with supporting legs	150 kg (331 lbs)

Colours: White Volcano grey

Material: Seat: polypropylene, Arm supports: polypropylene and TPE, Legs: aluminium and stainless steel, Ferrules: TPE





Etac – setting standards for more than 40 years

Scandinavian design and innovation are the foundation. Add to that quality, simplicity, comfort and functionality. For more than forty years Etac has developed products and services that create possibilities. Our products are developed in close cooperation with users and industrial designers. Welcome to the wonderful world of Etac!

Etac AB
Kista Science Tower
SE-164 51 Kista, Sweden
Tel +46 371 – 58 73 30
Fax +46 371 – 58 73 90
info@etac.se www.etac.com

etac[®]
Creating Possibilities